

A large, dark blue rectangular area containing the title text in white, sans-serif font. The background of the entire page is a close-up, high-angle photograph of a cryopreservation freezer's control panel, showing a digital display with '-80' and various buttons and indicators.

## CRYOPRESERVATION AND ULTRA LOW TEMPERATURE SOLUTIONS

PHCbi's VIP Ultra Low Temperature Freezers represent the industry's most complete combination of refrigeration technology, control, alarm, monitoring and accessibility for reliable sample preservation.

### Model Ranges:

MDF Cryogenic Freezers | CBS Cryogenic preservation  
TwinGuard ULT Freezers | FrostLess Freezers  
VIP ECO SMART Freezers | VIP ECO Freezers  
PRO ECO ULT Freezers




The complete line of PHC Europe B.V. products includes an array of laboratory equipment with advanced technology, controls, construction and performance attributes in the industry. Today we apply sophisticated refrigeration compressor design and state-of-the-art electronics to ultra-low and cryogenic freezers marketed to life science, pharmaceutical, biotechnology, healthcare and industrial laboratory markets.

PHCbi core technologies, patents, and intellectual properties are represented in every product line.

Core technologies apply to critical components and processes such as compressors, microprocessor electronics, and patented VIP PLUS vacuum insulation panels. These are engineered to exact specifications for important applications in life science, pharmaceutical, biotechnology, clinical, and industrial laboratories. As a result, PHCbi products operate with dependability, safety, energy efficiency, and ergonomic sensitivity.

# Your innovative partner in **ULTRA LOW TEMPERATURE** solutions

## Table of contents



|   |    |
|---|----|
| Commitment to quality                       | 4  |
| Cryopreservation                            | 6  |
| <b>MDF</b> Cryopreservation -150°C freezers | 8  |
| <b>CBS</b> Cryogenic preservation           | 12 |
| Ultra Low Temperature Freezers              | 18 |
| <b>TwinGuard</b> ULT Freezers (-86°C)       | 22 |
| <b>FrostLess</b> Freezers (-86°C)           | 28 |
| <b>VIP ECO SMART</b> Freezers (-86°C)       | 30 |
| <b>VIP ECO</b> Freezers (-86°C)             | 32 |
| <b>PRO ECO ULT</b> Freezers (-86°C)         | 36 |
| <b>Hybrid</b> technology                    | 38 |
| Cryogenic Vials                             | 42 |
| Racks                                       | 44 |
| Validation                                  | 48 |

## Committed to Quality

We at PHC Europe B.V. are committed to providing our customers with first class biomedical and pharmacy automation products and supporting services. Our goal is to maintain our market-leading reputation for excellent standards and for consistently living up to our customers' expectations.

Therefore, we are very proud that PHC Europe B.V.'s quality management system is certified (by TÜV Nederland) to **ISO9001:2015**

We hold the certification that covers the sales, distribution, service and validation of biomedical and pharmacy-automation equipment covering our sites in Etten-Leur (Netherlands), Avon (France) and Loughborough (UK).



**Life Science  
Innovator  
Since 1966**

Quality - It's a word that all are familiar with and one which many manufacturers claim to have or strive for. In the view of PHC, quality is a term that is ultimately defined by the customer. When quality becomes a customer-driven concept, quality includes meeting or exceeding our customers' needs or expectations. We focus on total quality which includes advanced processes and the culture of our organization. The result of our total quality initiatives involve many steps to provide our customers with superior value.

### 10 steps to guarantee superior value

#### Understanding Customer Demands

Our vision of total quality involves many face-to-face visits to customers to directly hear what they have to say about using the equipment we have supplied.

#### Creating a New Product Concept

Developing a concept for a new product is very similar to the "basic research" processes in life sciences. Technical staff and engineers develop various basic and innovative technologies to realize the concept for a new product.

#### Creating New Product Designs

As soon as engineering personnel begin developing a technological element, the design staff sets out to work on images.

#### Checking Local Legislation

We manufacture products for use in approximately 110 countries and regions around the world. Obviously, different laws and regulations apply, so we are always working to ensure our products conform to the laws and regulations of the individual locations

#### Design Review

Beyond determining specifications, evaluation criteria and achievement levels affecting product quality such as reliability, durability and safety standards must be achieved in mass-production models.

Founded in 1990 as subsidiary of the PHC Holdings Corporation, it is our mission to become a leading, trusted brand for sustainable healthcare and biomedical product solutions, which support the work of our customers to improve the health and well-being of people around the world.

For more than 30 years now, we respond to the needs of our pharmaceutical, biotechnology, hospital/clinical and industrial customers, offering an unique perspective on scientific research in general. As a result we play a critical role in product development for worldwide applications and have established a reputation as a manufacturer of high-quality and innovative medical and laboratory equipment.

Long lasting relationships have been built with leading pharmaceutical, healthcare and biotechnology companies as well as with major academic and research institutes in Europe. PHC Europe B.V. has set

the standard in many aspects. V.I.P. PLUS panels, Cool Safe compressors, Active Background Contamination Control and the world's first -152°C ULT freezer. Where PHC Europe B.V. took the initiative, the others followed. This made us a very important player in both the ultra-low temperature and the CO<sub>2</sub> market.

PHC Europe BV, part of the PHC Group, Biomedical Division, manages sales, marketing, logistics and technical service of PHCbi laboratory products throughout Europe, Africa, Middle East and parts of South America. Headquartered in the Netherlands with sales and service organizations in the UK, France and the Netherlands.

In the Dutch warehouse, about 1.000 units are waiting to be delivered directly from stock. Within a couple of days, spare parts can be delivered in every part of Europe. That's one of the strengths of PHC's European sales organization.



### Quality Assurance

Under PHC standards, quality actions required in new product development must clear three hurdles: AQ-0 approval for design completion, AQ-1 decision for shifting to mass production and AQ-2 decision for shipping mass-produced units.

### Mass Production

The Gunma factory was established in 1959 as the Tokyo manufacturing site for Sanyo Electric Co., Ltd. It is a core facility housing the Product Technology Development & Design Department, Quality Assurance Department and Production Department.

### Transportation packaging design

Our product packaging is designed to fulfill various distribution challenges around the world. Using past accumulated lessons learned we focus on reducing and recycling materials to support environmental initiatives.

### Installation

Product installations are often performed carefully by regional suppliers who are trained specialists. The path from facility entrance to the laboratory is measured and examined in advance to ensure precise, efficient installation. In some regions "white glove" service is also available to aid in installation and set-up.

### Sales and After-Sales Service

We always strive to learn more about our products. Before new products go into mass production at the factory, our sales personnel receives extensive knowledge about the products through intensive sales training

# Step into the world of Cryopreservation

Cryopreservation refers to the storage of a living organism, cell or tissue at ultra-low temperatures such that it can be restored to the same viable state as before it was frozen. Storage for an indefinite amount of time requires samples to be maintained below the glass transition temperature of aqueous solutions, approximately  $-130^{\circ}\text{C}$ , the temperature at which frozen water no longer sublimates and recrystallizes. Therefore  $-150^{\circ}\text{C}$  mechanical freezers or liquid nitrogen storage tanks are required for long-term preservation.

For years cryopreservation in liquid nitrogen at  $-196^{\circ}\text{C}$  has been the standard for long term preservation.

Over the years it became clear that cross contamination was a real risk with this kind of preservation. Not only can samples be contaminated by the nitrogen itself, but also by other samples. Evidence of cross contamination has been documented in the past, including a well known example of Hepatitis B transmission.

For this reason, storage in vapour phase became more popular. Valuable samples were no longer stored in liquid phase nitrogen to eliminate the risk of cross contamination.

There are also disadvantages of storage in vapour phase though:

1. Loss of storage capacity; The bottom third of the vessel cannot be used as there is liquid phase nitrogen there.
2. The temperature provided by vapour phase storage is highly dependent on the level of liquid nitrogen in the bottom of the vessel.
3. Large vertical temperature gradients can put valuable samples at greater risk during long term storage.

Dry storage systems, like  $-150^{\circ}\text{C}$  mechanical freezers and patented Isothermal storage systems, have been developed to offer the highest levels of protection for valuable samples, without the risks of cross contamination and vertical temperature gradients.

## Select your Cryopreservation solution

|  | V-Series<br>$-190^{\circ}\text{C}$ | Cryogenic Freezers<br>$-150^{\circ}\text{C} / -152^{\circ}\text{C}$ |
|--|------------------------------------|---|
| Storage in ULT-Freezer   | X                                  |   |
| Storage in dry vapour of $\text{LN}_2$                               | X                                  |   |
| Storage in vapour phase of $\text{LN}_2$                             |                                    |   |
| No risk of cross contamination through $\text{LN}_2$                 | X                                  | X   |
| Small temperature gradient risk                                      | X                                  | X   |
| Assured storage below critical temperature of $-130^{\circ}\text{C}$ | X                                  | X   |
| Autofill and monitoring required                                     | X                                  | not applicable  |
| Storage of more than 6000 2ml vials                                  | X                                  | X   |
| Ease of handling racks   | X                                  | X   |
| Usage of $\text{LN}_2$   | moderate                           | none  |
| No special requirements needed in storage room                       |                                    | X   |

# Cryogenic ULT -150°C freezers: dry storage

PHCbi has established a worldwide reputation as a manufacturer of high quality medical equipment over the past forty years. During this time we have become the leader in the European ultra-low temperature market. In this field PHCbi has set the standard in a number of ways including the introduction of VIP PLUS panels, application specific compressors and the world's first -152°C ULT freezer. Where PHCbi took the initiative, the others followed. Together with the high standard of service we deliver it makes us a major player in the biomedical market.

Due to the extended operation of compressors within ultra-low temperature freezers, lubricant oil is essential to reduce wear, prevent abrasions and seizure and therefore maintain the highest levels of freezer performance. However, when lubricant oil circulates within the refrigeration circuit, it can cause piping to become clogged, resulting in compressor damage.

Incorporating a high-efficiency oil separator, **Cryogenic** -150°C freezers effectively separate lubricant oil from refrigerant, increasing the durability of the compressors and offering a reliable ultra-low temperature environment. Highly efficient compressors have been specially developed and incorporated our **Cryogenic** -150°C freezers.



## -150°C Cryogenic Freezers are already in use to store:

- Umbilical cord blood as a source of hematopoietic stem and progenitor cells.
- Stem cells for autologous transplants in patients who have undergone high dose chemotherapy.
- Adipose tissue, epithelial cells and bone marrow for stem cell therapy.
- Blood products for immunology analysis.
- Mesenchymal stromal cells for regenerative medicine and tissue engineering.
- Cancerous tissue samples.
- Semen for artificial insemination – used for breeding guide dogs and race horses.
- Oocytes and embryos for IVF.
- Ovarian tissue for preserved reproductive function in women undergoing treatments.
- Plant seeds/shoots for breeding.
- Components for shrink fitting of parts in industrial manufacturing.
- Temperature testing of parts in the automotive and aerospace industry.

# MDF-C2156VAN Cryogenic Freezer

**Cryogenic Freezers** are well-known for maintaining uniform temperatures at  $-150^{\circ}\text{C}$  for the reliable, long-term preservation of cells and tissues. With thin vacuum insulation panel (VIP PLUS) walls, the MDF-C2156VAN **Cryogenic Freezer** can achieve more storage capacity than a conventionally insulated freezer without increasing footprint, while also maintaining superior temperature uniformity.

## VIP PLUS INSULATION



PHCbi's patented VIP PLUS technology has resulted in a revolutionary vacuum insulation cabinet construction with improved thermal properties for superior temperature performance.

## CRYOGENIC FREEZERS PROVIDE SAMPLE STABILITY

A uniformity of  $\pm 5^{\circ}\text{C}$  in our mechanically refrigerated **Cryogenic Freezers** is far superior to the top-to-bottom temperature uniformity provided by liquid nitrogen vapour phase storage, without the concern of cross-contamination often associated with liquid nitrogen (liquid phase storage).

When you need long-term preservation at  $-150^{\circ}\text{C}$ , put your trust in PHC Cryogenic freezers.

- 230V / 50Hz connection supply for flexible use in the laboratory, without additional installations
- Low heat emission. Allows multiple devices to be placed in one room
- Cool-Safe compressors increase the reliability of long-term preservation
- High performance cooling with low noise level designed by PHCbi
- Alarm and safety functions
- Remote alarm contact (NO / NC)
- Multiple access ports allow independent probe-scan be introduced
- Emergency cooling with liquid  $\text{N}_2$  already standard on some models



Model: MDF-C2156VAN-PE



## ULTRA-LOW -150°C FREEZER WITH VIP PLUS INSULATION

**Cryogenic** Freezers with VIP PLUS vacuum insulation provide up to 30% more storage capacity than a conventionally insulated freezer, without increasing the footprint. A glass fibre core provides advanced thermal properties. This results in a large capacity -150°C freezer for storage of up to 150 world standard 2" boxes with a minimal footprint.

## CONTROL PANEL WITH GRAPHIC LCD DISPLAY


All alarm functions, self-diagnostic notifications and a graphical display of temperature performance over time are available in the specially designed LCD control panel. The blue display provides a clear view of the temperature and gives a notification in the case of abnormalities in temperature, ambient temperature, power supply etc.

## APPLICATION SPECIFIC COMPRESSORS

The MDF-C2156VAN is equipped with compressors that are specifically designed for ultra-low temperature applications. These compressors achieve a 10% reduction in energy consumption and the aerodynamically designed and placed components in the refrigeration compartment provide superior airflow, significantly reducing the stress to the freezer and contributing to excellent durability.




### Liquid Nitrogen freezer (Traditional vapour phase)



### MDF-C2156VAN

mechanically refrigerated **Cryogenic** Freezer



Comparison of temperature distribution in a liquid nitrogen freezer (vapour phase) and our MDF-C2156VAN mechanically refrigerated **Cryogenic** Freezer. The graph shows temperatures at different locations within the chamber. This data demonstrates that 100% of the MDF-C2156VAN storage space maintains uniform storage temperatures safely below -130°C, while temperature in the LN<sub>2</sub> vapour system is dependent on storage location.

# MDF-1156 Cryogenic Freezers

The Cryogenic MDF-1156-PE Ultra Low Temperature Freezer is a mechanically refrigerated cryogenic freezer without harmful CFCs. This Ultra Low Temperature Freezer achieves stable long-term preservation of cells and tissues. Highly efficient compressors have been specially developed and incorporated in the freezing unit. With a powerful low noise design afforded by traditional ultra-low temperature technology, this freezer delivers durable, stable cooling.


## MICROPROCESSOR CONTROL WITH LED DIGITAL DISPLAY

Accurate temperature setting and a stable performance are provided by a microprocessor controller with an easy to use LED digital display and flat key data entry. The MDF-1156 maintains an ultra-low temperature of  $-152^{\circ}\text{C}$ , far lower than the recrystallization point for pure water ( $-130^{\circ}\text{C}$ ). This low temperature provides the ideal preservation environment for long-term storage.

## ALARM AND SAFETY DEVICES ENSURE HIGH LEVELS OF PROTECTION

A high temperature alarm alerts users when the temperature of the freezer deviates  $10^{\circ}\text{C}$  or  $15^{\circ}\text{C}$  above the set temperature and a power failure alarm is activated if power to the freezer is lost. As well as lamp and buzzer functions, a remote alarm contact can be set up to alert users to these abnormalities outside of working hours.

## RECRYSTALLIZATION MECHANISM



Model: MDF-1156-PE




| Cryogenic Freezers          |          |                  |                   |
|-----------------------------|----------|------------------|-------------------|
| Model Number                |          | MDF-1156-PE      | MDF-C2156VAN-PE   |
| Temperature control range   | °C       | -150             | -150              |
| External dimensions (WxDxH) | mm       | 1400 x 800 x 945 | 1730 x 765 x 1010 |
| Internal dimensions (WxDxH) | mm       | 500 x 450 x 572  | 760 x 495 x 615   |
| Volume                      | litres   | 128              | 231               |
| Capacity                    | 2" boxes | 81               | 150               |
| Net weight (approx)         | kg       | 265              | 318               |

\* Complies with Art. 11, Annex III of F-Gas Regulation (EU) No 517/2014. Contains fluorinated greenhouse gases.

## Optional accessories

| OPTIONS                       |  |             |                      |
|-------------------------------|--|-------------|----------------------|
| Cryogenic Freezers            |  |             |                      |
| Model Number                  |  | MDF-1156-PE | MDF-C2156VAN-PE      |
| Liquid N <sub>2</sub> back-up |  | CVK-ATN2-PW | Supplied as standard |
| Temperature recorders         |  |             |                      |
| - Continuous strip type       |  | MTR-155H-PW | MTR-155H-PW          |
| - Chart paper                 |  | RP-155-PW   | RP-155-PW            |
| - Ink pen                     |  | DF-38FP-PW  | DF-38FP-PW           |
| - Recorder housing            |  | -           | MDF-S30150-PW        |



# CBS Isothermal Freezers: No liquid nitrogen contact

Liquid nitrogen is, like all liquids, a possible transporter of contaminants. Contamination between samples, leakage of liquid nitrogen into sample vials during storage and contact of liquid nitrogen with skin are just some of the risks that can be eliminated with an Isothermal liquid nitrogen dry storage freezer from Custom BioGenic Systems (CBS).

CBS Isothermal freezers feature a patented liquid nitrogen jacket to provide uniform storage temperatures in the  $-190^{\circ}\text{C}$  range, free from liquid nitrogen contact.

In the year 2000 Custom Biogenic Systems patented the first  $-190^{\circ}\text{C}$  dry storage system. This new storage system uses the same vacuum insulated vessel as conventional  $\text{LN}_2$  freezers, but instead of the liquid nitrogen being in the sample storage space of the freezer, it is in a patented liquid nitrogen jacket in the wall of the freezer. This liquid nitrogen jacket produces a unique vapour circulation in the freezer. This vapour circulation, in combination with the convection from the freezer walls, produces a very stable temperature in the  $-190^{\circ}\text{C}$  range without the large vertical temperature

gradients found in traditional vapour phase  $\text{LN}_2$  storage tanks. The minimal temperature gradient and lack of liquid nitrogen in the bottom of the vessel means that the full internal volume of each freezer is available for storing samples at liquid nitrogen temperatures, making the most efficient use of capacity.

## Medical Device Regulation



All Isothermal models are certified to Medical Device Regulation (EU) 2017/745

### THE ISOTHERMAL CONCEPT

The sample storage area is cooled by a liquid nitrogen jacket surrounding the stainless steel interior, and by nitrogen vapour entering the freezer from the jacket via directional vents. This patented technology provides exceptional temperature uniformity in the  $-190^{\circ}\text{C}$  range, allowing the full freezer capacity to be used with confidence. The circulation of vapour within the freezer also results in less cold air loss during lid opening and improved visibility. This allows full-width lids to be used providing quick, unrestricted access to sample racks.

### NO LIQUID NITROGEN CONTACT

With no liquid nitrogen in the storage area, samples can be stored safely in the  $-190^{\circ}\text{C}$  range without the risk of cross-contamination through liquid nitrogen. The Isothermal design also provides added user safety with no splashing or contact of liquid nitrogen when removing racks and samples from the freezer.



The innovative design of CBS Isothermal freezers eliminates a number of major risks of traditional LN<sub>2</sub> storage, including:

- **Cross contamination:** Studies have shown that viral, bacterial and fungal pathogens can survive after suspension in liquid nitrogen. Infected samples can cross contaminate other samples in the same liquid nitrogen tank.
- **Loss of samples:** Storage directly in liquid nitrogen can make vials shrink. This may cause liquid nitrogen to seep into the vials, which on rewarming, expand and subsequently explode as nitrogen vapourizes inside the vials.
- **Health and safety:** Traditional storage in liquid nitrogen exposes users to direct contact with LN<sub>2</sub> which can lead to cold skin burns.

CBS Isothermal freezers feature a patented liquid nitrogen jacket to provide uniform storage temperatures in the -190°C range, free from liquid nitrogen contact.

## AUTOMATIC OPERATION

Isothermal freezers feature the series 2301 auto-fill and monitor system, which controls the automatic filling of the liquid nitrogen jacket and provides the user with an easy to read overview of the freezer temperature and status.

## SAMPLE SECURITY

A comprehensive alarm system with remote alarm contact constantly monitors all aspects of the freezer's operation. Samples are also protected by lid and control panel locks. The freezer can be monitored by a central BMS or monitoring system.

## IMPROVED USABILITY

The 2301 SmartTouch Controller user interface, designed specifically for added peace of mind, features dual temperature display, autofill, fill timer, exportable reports, Modbus enabled, multiple alarms and is monitoring system ready.

# Rack configurations

Organizing the space in your freezer can help you to make your work more efficient. Good sample management can result in both cost and time savings. By choosing the right racks, not only will efficiency be increased, but the risk of sample degradation due to exposure to ambient temperatures will be significantly reduced. This is a fact, whether you store stem cells, cord blood, bone marrow or any other type of cell or tissue sample.

Whether your storage unit a chest freezer or even a liquid Nitrogen freezer, an organized freezer will provide you with:

- Time efficiency because you can locate, retrieve and replace your samples easily and quickly.
- Cost efficiency because organized samples and cell lines might reduce the number of freezers.
- Safety because your samples are better protected.



# Straw storage inventory system

Rotating carousels provide access to lower storage levels and maximise storage space.

The upper level provides a convenient workspace where canisters can be parked (shown in green) and samples inspected at -190°C.

Use the retrieval tool to place and retrieve canisters quickly and easily.



Isothermal -190°C Dry Storage Freezers

| Model Number                     |        | V-1500AB-ST      | V-3000AB-ST       | V-3000ABEH-ST     | V-5000AB-ST        | V-5000ABEH-ST      |
|----------------------------------|--------|------------------|-------------------|-------------------|--------------------|--------------------|
| Liquid nitrogen capacity         | litres | 30               | 70                | 89                | 93                 | 140                |
| <b>Dimensions</b>                |        |                  |                   |                   |                    |                    |
| External dimensions (W x D x H)  | mm     | 660 x 939 x 1143 | 939 x 1219 x 1206 | 939 x 1219 x 1473 | 1219 x 1371 x 1320 | 1219 x 1371 x 1473 |
| Usable interior height           | mm     | 736              | 736               | 940               | 736                | 864                |
| Usable interior diameter         | mm     | 534              | 787               | 787               | 1016               | 1016               |
| Weight empty                     | kg     | 148              | 272               | 295               | 425                | 453                |
| Weight full                      | kg     | 174              | 327               | 367               | 500                | 566                |
| <b>Maximum capacity</b>          |        |                  |                   |                   |                    |                    |
| Max. vial capacity (2ml)**       | qty    | 9100             | 22100             | 25500             | 40300              | 46500              |
| Max. blood bag capacity (50ml)** | qty    | 434              | 1120              | 1280              | 1936               | 2208               |

\*\* Capacity is subject to rack type

# Optional accessories



Isothermals, LN<sub>2</sub> freezers & cryosystems

|                                 |   |                               |  |
|---------------------------------|---|-------------------------------|--|
| Cryokit Glove 400               | See below   | LN <sub>2</sub> Level stick   | <ul style="list-style-type: none"> <li>• 1/2 Centimeter and 1/4 inch increments.</li> <li>• Can withstand temperatures up to -190°C.</li> <li>• Measures up to 36" [92cm].</li> </ul>  |
| "T" Valve                       | Solid brass cryogenic shut-off valve (rated for temperatures from -196° C to 74° C). <ul style="list-style-type: none"> <li>• 2 Male 1/2" NPT brass fittings.</li> <li>• 1 Female 1/2" NPT stainless steel Flared fitting.</li> </ul> | LN <sub>2</sub> Transfer hose | <ul style="list-style-type: none"> <li>• Flexible stainless steel construction.</li> <li>• 1/2" NPT flared fitting on both ends. (3/8" I.D.)</li> <li>• 4', 6' lengths are available (custom lengths are available upon request).</li> </ul> |
| "Y" Valve                       | <ul style="list-style-type: none"> <li>• 2 Male 1/2" NPT brass fittings.</li> <li>• 1 Female 1/2" NPT stainless steel flared fitting.</li> <li>• Overall length approximately 6".</li> </ul>  | Roller base                   | Cryosystem roller bases. Reinforced fiberglass construction with casters. 5 configurations available   |
| LN <sub>2</sub> Phase separator | Designed to minimize hazardous splashing and vapourization, phase separators are available to use when transferring liquids into various open containers.   | Canes                         | <ul style="list-style-type: none"> <li>• 2ml cane.</li> <li>• 5 and 6 place vial canes.</li> <li>• Lightweight aluminum construction.</li> <li>• Overall length approximately 11.5" (6 place cane).</li> </ul>                               |
| Cardboard sleeves               | <ul style="list-style-type: none"> <li>• 5 and 6 place sleeves for standard 2ml cane.</li> <li>• Cardboard construction.</li> </ul>   | Hose covers                   | <ul style="list-style-type: none"> <li>• Water resistant fabric 4', 6' lengths are available.</li> </ul>   |

## Cryokit Glove 400

CRYOKIT GLOVE is a fully waterproof glove to protect the hand and the arm in operations that involve the use of Liquid Nitrogen (or other cryogenic liquids and gases).

Made with specially designed composite materials for this application, it keeps the hands warm and dry even for long periods of use. The CRYOKIT GLOVES are proven effective to protect the hands and arms, by contact with Liquid Nitrogen (-195.82 °C) for one minute and maintain their flexibility (RCT Test method). The multilayered construction allows a high level of thermal protection, flexibility and dexterity at -200 °C.

Indicated for all work in the presence of Liquid Nitrogen or other liquids and gases, to prevent contact of cold and burns in the event of leakage of the cryogenic liquid.

The glove is made entirely with specially coupled composite materials. On the inside there is a integral membrane which allows perspiration but does not allow liquids to penetrate inside.

The cold protection is ensured by an inner lining made of soft fleece, which also coupled with polyester.



CE certificate in 3rd Cat.

# IntelliRate™ i67C

Perfectly sized to efficiently scale precision rate freezing as your needs grow.

As cell-based research and bioproduction volumes increase, you need a repeatable freezing process that can efficiently improve post-thaw viability rates and cryopreservation quality at scale. CBS designed the new IntelliRate i67C to increase single-batch, liquid nitrogen (LN2) controlled rate freezing (CRF) capacities to meet your needs as development and manufacturing volumes increase.

## Increase Rate Freezing Batch Capacity and Quality Cost Effectively

Expand controlled rate freezing throughput with significantly larger batches, without adding significant cost. With 40% more capacity than any other benchtop CRF, not only can the IntelliRate i67C increase production levels, but also ensure more process consistency and time savings by running thousands of samples in a single run.

## Use Workspace and Staff More Efficiently

With its compact footprint accommodating larger capacity, the IntelliRate i67C frees up valuable workspace by using a single unit to freeze volumes, where normally multiple standard size freezers would be required. Fewer runs also simplify programming, setup, LN2 monitoring and power requirements.



## Broaden Your Rate Programming Capabilities and Options

Use the largest touchscreen controller in its class to execute a wide range of preset and variable programming capabilities, including setup of repeatable, 21 CFR Part 11 compliant freeze profiles. With the IntelliRate i67C, a range of options are also available for enhanced multi-probe monitoring, racking and sample management.

| IntelliRate™ i67c              |              |  |
|--------------------------------|--------------|--|
| <b>Dimensions</b>              |              |  |
| External dimensions (WxDxH)    | mm           | 495 x 781 x 787  |
| Internal dimensions (WxDxH)    | mm           | 406 x 406 x 406  |
| Volume                         | litres       | 67   |
| Net weight                     | kg           | 49,9   |
| Capacity                       | 50mL bags    | 104  |
|                                | 2mL vials    | 2166   |
|                                | cassettes    | 100  |
|                                | 1-liter bags | 26   |
|                                | canes        | 272  |
| <b>Electrical</b>              |              |  |
| Power Supply                   | V            | 120/230  |
|                                | Hz           | 50/60  |
|                                | Amp          | 10   |
| <b>Programming and Control</b> |              |  |
| Preset programs                |              | 6 with unlimited variations  |
|                                |              | Repeatable temperature rate changes within cycles  |
| Programmable temperature range | °C           | -180°C to +50°C  |
| Freeze rates                   | °C           | 0.01 to 99.9°C per minute  |
|                                |              | 21 CFR Part 11 compliant reporting protocols   |
| <b>Options and Accessories</b> |              |  |
|                                |              | Multi-probe version provides temperature monitoring at six additional customer-specified positions                             |
|                                |              | Innovative racking systems: choose how you want to freeze, retrieve and move samples to long-term storage or into clinical use |
|                                |              | Broad array of LN2 management, monitoring, sample handling and freezer accessories   |



# 2101 Controlled rate freezer

The 2101 controlled rate freezer meets the highest standards for the programmed freezing of biological samples. Computer controlled temperatures ensure your samples are frozen at precisely the same rate during each run. Freeze protocols can be controlled by either the chamber or the sample temperature. The freezing rate can be programmed at the optimum rate for each individual sample.

The 2101 controlled rate freezer comes equipped with a dedicated laptop and 2100 programming software as standard to provide the highest levels of flexibility and user convenience.

**Benefits include:**

- Unlimited programming capability.
- Multi-colour graph for sample, chamber and program temperature.
- Sample or chamber temperature set-point control.
- Continuous control status indication.
- Programs and freeze data saved to hard drive or disc.
- Password protected software.
- Selectable password levels.
- Searchable database for freeze run history.

- Freeze run graphs and data available via any standard computer printer.
- 30 data field available for each freeze run.
- Continuous digital and graphical display of time and temperature during operation.
- Audible and visual indicators for:
  - End of each freeze run.
  - Temperature probes.
  - Freeze run tracking.
- On screen help.

**Standard system Includes:**

- Laptop computer with 2100 software.
- Freezing chamber.
- Choice of sample rack.
- 1.2 m LN<sub>2</sub> transfer hose.

**Options:**

- Choice of sample racks for vials, tubes, straws, bags & canes.
- Temperature probes for different sample types.
- Roller cart.



Model: 2101 controlled rate freezer

| 2101 Controlled Rate Freezer    |               |                 |                                   |        |   |
|---------------------------------|---------------|-----------------|-----------------------------------|--------|---|
| Dimensions                      |               |                 | Control                           |        |   |
| External Dimensions (W x D x H) | mm            | 484 x 648 x 770 | Controller                        |        | Laptop controller (included)                                  |
| Internal Dimensions (W x D x H) | mm            | 356 x 243 x 349 | Operating System                  |        | Windows based   |
| Volume                          | litres        | 28              | Temperature sensor                |        | type T thermocouple (2 as standard, up to 8 with multi-probe) |
| Net Weight                      | kg            | 34,7            | <b>Construction</b>               |        |   |
| Capacity                        | 1.2-2ml vials | 650             | Exterior Material                 |        | Powder coated stainless steel                                 |
|                                 | 4-5ml vials   | 390             | Interior Material                 |        | Stainless steel   |
|                                 | bag canisters | 10 - 20         | Max. load - total                 | kg     | 20  |
|                                 | canes         | 130             | <b>Electrical and Noise Level</b> |        |   |
| Programmable Temperature Range  | °C            | +50 to -180     | Power Supply                      | V      | 230   |
| Programmable Cooling Rate Range | °C / min      | 0.01 to 99.9    | Frequency                         | Hz     | 50  |
|                                 |               |                 | Noise Level                       | dB [A] | <66   |

# Ultra Low Temperature Freezers

In response to the needs of leading pharmaceutical, biotechnology, hospital/clinical and industrial customers, PHC Europe B.V. offer a unique perspective on Ultra Low Temperature preservation. As a Member of the PHC Group, Biomedical Division, PHC Europe B.V. plays a critical role in product development for worldwide applications.

As a result, PHCbi products incorporate the latest applied theories and ideas, refrigeration compressors, electronic components, energy efficiency, robotic manufacturing and economies-of-scale that directly benefit our customers. Furthermore, PHCbi products are extensively tested to meet the toughest quality standards in the world; our own! Through the vigorous application of our unique Vertical Component Integration™, PHCbi products offer tangible benefits - from performance and reliability to ergonomics and convenience, that no other manufacturer can provide.

## PHCbi DESIGN & ULTIMATE RELIABILITY

VIP ultra low temperature freezers offer advanced cabinet design, reliable refrigeration systems and easy-to-use controllers making them ideal for the long-term secure storage of valuable samples. Every component is carefully selected and matched for optimum operation under demanding laboratory conditions, while the internal layout of the refrigeration system is meticulously designed for maximum heat removal, reducing stress on the system and therefore providing the highest levels of reliability and durability.

## SUPERIOR PERFORMANCE

All PHCbi freezers are designed to provide the highest quality construction with superior performance. Key features such as strategically placed evaporator coils, V.I.P. PLUS panels and insulated inner doors contribute to the unrivalled temperature uniformity and stability of VIP freezers, allowing the freezers to conform to the strictest standards and validation protocols. Quieter operation is achieved through condenser fan blade design, noise reduction insulation, anti-vibration systems and internal compressor noise reduction.

## EFFICIENT COOLING

Cascade refrigeration systems within the VIP upright freezers provide efficient cooling with optimized heat exchange pathways and increased cooling capacity for reliable sample protection and cost effective operation at ultra low temperatures.

## ENHANCED USE & INTELLIGENT SECURITY

Our freezers are managed and monitored by an integrated microprocessor controller with a comprehensive alarm system and diagnostic functions. A Status Alert feature constantly monitors ambient and system conditions and notifies the user of any abnormalities before a problem occurs. The MDF-DU901VHL is updated with a touch screen that allows full user control, even with gloved hands, and a USB port for convenient transfer of logged data to a PC.

## Medical Device Directive



The MDF-DU502VH-PE, MDF-DU702VH-PE and MDF-DU300H-PE are certified as a Class IIa Medical Device (93/42/EEC and 2007/47/EC).

Applicable countries: Austria, Belgium, Cyprus, Denmark, Finland, France, Germany, Ireland, Italy, Liechtenstein, Luxembourg, Malta, the Netherlands, Spain, Switzerland and the United Kingdom only

For laboratory use

Applicable countries: EEA countries, Switzerland and Turkey


## VIP PLUS INSULATION



PHCbi's patented VIP PLUS technology has resulted in a revolutionary vacuum insulation cabinet construction with improved thermal properties for superior temperature performance.


## SUPERIOR FOOTPRINT

PHCbi ultra low temperature freezers with space-saving VIP insulation offer outstanding energy efficiency, whilst delivering exceptional cooling performance and durability for storing valuable research and clinical samples.



## INNOVATIVE DESIGN


PHC was the first company to introduce vacuum insulation panels to ultra low temperature and cryogenic freezers. The PHCbi patented V.I.P. PLUS vacuum insulation panel thin-wall composite is a high-efficiency design that yields more interior storage volume in a conventional freezer footprint. The PHCbi VIP Freezer range typically provide 30% more storage capacity for a given floor area saving valuable laboratory space.



## OPTIMUM UNIFORMITY

Uneven interior temperatures can lead to a loss in sample integrity. PHCbi freezers with uniform, stable temperatures and quick recovery times provide the best protection for your samples, ensuring reliable preservation while guarding against degradation.

Surpasses the customer preference of  $\pm 5^{\circ}\text{C}^*$



\* Based on internal validation data tested at  $-80^{\circ}\text{C}$  setpoint, in an empty chamber with  $23^{\circ}\text{C}$  ambient temperature.

\* The data may vary depending on the use, circumstances and optional accessories. Validation documents can be provided for each serial number for an additional fee.

## EXCELLENT RELIABILITY AND UNIFORMITY WITHIN AN OPTIMAL FOOTPRINT.

The refrigeration systems within PHCbi VIP freezers are designed specifically for demanding ultra low temperature applications to give proven durability. Manufactured with space-saving VIP PLUS vacuum insulation panels, they are ideally suited for use in laboratories and hospitals for long-term preservation of samples. Whatever your preservation needs are, PHCbi provides the right equipment for your application.

- 1 Multiple access ports permit insertion of independent probes, instrumentation or liquid CO<sub>2</sub> back-up injectors.
- 2 Universal keyed door lock offers added security.
- 3 A vacuum release port (available on the MDF-DU502VH-PE, MDF-DU702VH-PE and MDF-DU901VHL-PE) allows smooth door opening when the door seal is tightened by negative pressure generated from temperature difference between chamber and ambient.
- 4 Insulated and gasketed inner doors seal inside to offer additional protection, improve uniformity.
- 5 EZlatch for smooth, one-handed operation and positive gasket seal. Provision for padlock.
- 6 Temperature recorder (optional) mounts easily in pre-engineered mounting space.
- 7 PHCbi designed compressors are specifically designed for ultra low temperature applications.
- 8 High impact, recessed casters and leveling feet simplify installation.
- 9 An integrated microprocessor controller with LCD Touch Screen to simplify all freezer functions.
- 10 Front access to washable, electrostatic condenser filter for routine condenser air filter cleaning.
- 11 Heated outer door gaskets and a 'hot line' circulating hot refrigerant gas around the door frame ensure minimal ice build-up.



## VIP SERIES MONITORING



### CONTROLLED FREEZER ACCESS

(available on the MDF-DU503VH-PE and MDF-DU703VH-PE)


Know who entered the freezer and when.  
Track how long the door was open.

**4 Ways to Lock**

- Facial recognition
- Keypad
- Proximity card (Optional)
- Unique individual key



electronic door lock




facial recognition

## VACUUM RELEASE PORT

A vacuum release port allows smooth door opening when the door seal is tightened by negative pressure generated from temperature difference between chamber and ambient.

## FLEXIBLE SHELF LAYOUT

Multiple shelf configurations in the upright models allow a variety of storage options. Organize your samples by using your existing inventory racks or select from the many different rack types we offer. PHCbi's racks are made of stainless steel or anodized aluminum. The aluminum racks are very light, yet sturdy and corrosion free.



## EZlatch



The EZ Latch door handle was developed based on human engineering. It is designed for secure door opening/closing with minimum effort as the name suggests and make access to stored sample easier.

# TwinGuard ULT Freezers

**TwinGuard** Ultra Low Temperature Freezers with Dual Cooling Technology offer the highest level of security for high-value samples. Alongside exceptional ease-of-use and data monitoring, the Dual Cooling System provides the highest level of protection through the use of two independent refrigeration systems. If one system unexpectedly fails, the other can maintain the freezer's temperature uniformly in the  $-70^{\circ}\text{C}$  range. Developed for use with conventional inventory racks and boxes, the **TwinGuard** Series is ideal for storage of sensitive, high-value samples.

**When sample security and peace of mind are of paramount importance, put your trust in TwinGuard ultra low freezers.**


- An optimal Dual Cooling System provides an unparalleled level of safety and added peace of mind through the use of two independent refrigeration systems.
- ECO mode overlaps refrigeration cycles to reduce energy consumption.
- Advanced space-saving VIP PLUS technology (compared to our conventional models).
- Graphical LCD display with data monitoring and data log exported by USB.
- Available in Upright and Chest models.
- Filterless design reduces the time for routine maintenance
- Two Vacuum release ports
- No icing on frame
- Alarm and safety functions

## SCIENTIFIC APPLICATIONS

- Temperature sensitive samples such as therapeutics and biospecimens.
- Samples needing to retain viability such as stem cells, engineered tissue, organs, vaccines, hydromas, cancer cells or fibroblasts.
- Longitudinal study samples.
- Important medical research samples.
- Valuable pharmaceutical products.
- Clinical trial samples.
- Pathogenic samples within high security laboratories.




Model: MDF-DC700VX-PE



Model: MDF-DU702VX-PE

## DUAL COOLING SYSTEM



Within *TwinGuard's* independent Dual Cooling System, efficient ultra-low cooling is achieved through two independent evaporator circuits surrounding the interior chamber.

## Medical Device Directive



The MDF-DU302VX-PE, MDF-DU502VX-PE, MDF-DU702VX-PE, MDF-DC500VX-PE and MDF-DC700VX-PE series are certified as a Class IIa Medical Device (93/42/EEC and 2007/47/EC).

Applicable countries: Austria, Belgium, Cyprus, Denmark, Finland, France, Germany, Ireland, Italy, Liechtenstein, Luxembourg, Malta, the Netherlands, Spain, Switzerland and the United Kingdom only


For laboratory use

Applicable countries: EEA countries, Switzerland and Turkey

### DUAL COOLING SYSTEM

Within *TwinGuard's* independent systems, efficient ultra-low cooling is achieved through two independent evaporator circuits surrounding the interior chamber.

### Two independent evaporator circuits



Dual Cooling System Upright freezers

### ENHANCED USE & INTELLIGENT SECURITY

The freezers are managed and monitored by an integrated microprocessor controller with a comprehensive alarm system and diagnostic functions. Status and control of parameters are accessible via an LCD information centre.

The EZlatch, on the upright models, makes access to stored samples even easier. A colour LCD touch panel allows full user control, even with gloved hands, while the USB port makes transferring logged data to a PC convenient.



- The Dual Cooling System offers the highest level of security through the use of two independent refrigeration systems. If one system unexpectedly fails the other can maintain the freezer at the  $-70^{\circ}\text{C}$  range.

In the case of unexpected failure of one of the cooling circuits, the other circuit will maintain the freezer continuously in the  $-70^{\circ}\text{C}$  range.

## MEETING YOUR FREEZER STORAGE NEEDS

An organized freezer will provide you with:

- Time savings - locate, retrieve and replace your samples easily and quickly.
- Cost savings - organized samples and cell lines can help to reduce the number of freezers.
- Added sample security and energy savings - samples are better protected and are less exposed to ambient temperatures as door opening times can be reduced when placing and retrieving samples, which also reduces energy use.

For an overview of the racks designed for the TwinGuard series see page 42-45.

## FILTERLESS DESIGN

The filterless construction of the freezers reduces routine maintenance time by eliminating the need for regular cleaning of filters.

## FLEXIBLE SHELF LAYOUT

Multiple shelf configurations in the upright models allow a variety of storage options. Organize your samples by using your existing inventory racks or select from the many different rack types we offer. PHCbi's racks are made of stainless steel or anodized aluminum. The aluminum racks are very light, yet sturdy and corrosion free.

- High quality racks – designed for safe working and easy access to samples.
- Affordable solutions – making freezer storage cost-effective as well as space-efficient.
- Large selection of products – additional rack types and boxes are available on request.


## SUPERIOR FOOTPRINT

PHCbi ultra-low temperature freezers with space-saving VIP PLUS insulation offer outstanding energy efficiency, whilst delivering exceptional cooling performance and durability for storing valuable research and clinical samples.

### PRESERVE SAMPLE INTEGRITY FOR BETTER END PRODUCTS


Uneven interior temperatures can lead to a loss in sample integrity. Freezers with uniform, stable temperatures and quick recovery times provide the best protection for your samples, ensuring reliable preservation while guarding against degradation.

FIGURE 2A - MDF-DU702VX; 9 POINT TEMPERATURE MAPPING



### -86°C ULTRA-LOW FREEZERS (TwinGuard)

FIGURE 2B







## INDUSTRY FIRST INTELLIGENT ECO MODE OPERATION

The **TwinGuard** series freezer can be set to Normal or ECO mode operation, depending on the requirements of the user. Although both refrigeration systems are completely independent, ECO mode establishes an overlapping cycle to significantly reduce energy consumption while maintaining optimum interior uniformity for protection of high value materials. Normal mode maintains the most repeatable, cycling wave form for the strictest of GMP applications.

### ECO MODE OPERATION IMAGE



### NORMAL MODE OPERATION IMAGE



TwinGuard Upright Freezers

| Model Number                | MDF-DU302VX-PE                    | MDF-DU502VX-PE                      | MDF-DU702VX-PE                      |
|-----------------------------|-----------------------------------|-------------------------------------|-------------------------------------|
| Temperature control range   |                                   | -50 --86                            |                                     |
| External dimensions (WxDxH) | 670 x 882 x 1840                  | 790 x 882 x 1993                    | 1030 x 882 x 1993                   |
| Internal dimensions (WxDxH) | 490 x 600 x 1230                  | 630 x 600 x 1400                    | 870 x 600 x 1400                    |
| Volume                      | 360                               | 528                                 | 729                                 |
| Capacity                    | 240                               | 384                                 | 576                                 |
| Power Consumption           | Normal Mode: 9,7 / Eco Mode: 8,3* | Normal Mode: 16,5 / Eco Mode: 15,3* | Normal Mode: 16,9 / Eco Mode: 15,7* |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load]

\* Complies with Art. 11, Annex III of F-Gas Regulation (EU) No 517/2014. Contains fluorinated greenhouse gases.



### Touch Panel Legend


1. **Present temperature display field:**  
The current chamber temperature is displayed.
2. **Set temperature value display field:**  
The set value of chamber temperature is displayed.  
Default setting: -80°C.
3. **Message display field:**  
Alarms, errors or messages are displayed when a fault occurs.
4. **Control display:**  
The present operation control mode is displayed.  
Normal control: "Normal" is displayed.  
Eco control: "ECO" is displayed
5. **Alarm display:**  
Normal condition: "Normal" is displayed.  
Alarm-activated, buzzer-delayed: "Alarm" is displayed.  
Alarm-activated, buzzer-sounding: "Warning" is displayed.
6. **Outer door (opening / closing display)**

### MICROPROCESSOR CONTROL WITH TOUCH SCREEN DISPLAY

The **TwinGuard** series is managed by an integrated microprocessor controller with LCD information center to simplify all freezer functions. Uniform ultra-low temperature is achieved through a combination of performance systems supervised by the controller complete with alarm, programming and diagnostic protocols.

The **TwinGuard** chest freezers feature a touchscreen display and a USB port which allows logged data to be easily transferred to a PC.

### DUAL COOLING SYSTEM IN THE VX CHEST FREEZERS





| TwinGuard Chest Freezers    |          |                                     |                                     |
|-----------------------------|----------|-------------------------------------|-------------------------------------|
| Model Number                |          | MDF-DC500VX-PE                      | MDF-DC700VX-PE                      |
| Temperature control range   | °C       | -50 ~-86                            |                                     |
| External dimensions (WxDxH) | mm       | 2010 x 845 x 1070                   | 2300 x 845 x 1070                   |
| Internal dimensions (WxDxH) | mm       | 1190 x 640 x 756                    | 1480 x 640 x 756                    |
| Volume                      | litres   | 575                                 | 715                                 |
| Capacity                    | 2" boxes | 416                                 | 520                                 |
| Power Consumption           | kWh/day  | Normal Mode: 15,5 / Eco Mode: 13,6* | Normal Mode: 16,3 / Eco Mode: 14,9* |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load]

\* Complies with Art. 11, Annex III of F-Gas Regulation (EU) No 517/2014. Contains fluorinated greenhouse gases.

## Optional accessories

| TwinGuard ULT Freezers         |          |                |                |                |                |                |
|--------------------------------|----------|----------------|----------------|----------------|----------------|----------------|
| Model number                   |          | MDF-DU302VX-PE | MDF-DU502VX-PE | MDF-DU702VX-PE | MDF-DC500VX-PE | MDF-DC700VX-PE |
| Liquid CO <sub>2</sub> back-up |          |                | MDF-UB7-PW     |                | MDF-UB7-PW     |                |
| Temperature recorders          |          |                |                |                |                |                |
| - Circular type                |          |                | MTR-G85C-PE    |                | MTR-G85C-PE    |                |
| - Chart paper                  |          |                | RP-G85-PW      |                | RP-G85-PW      |                |
| - Ink pen                      |          |                | PG-R-PW        |                | PG-R-PW        |                |
| - Continuous strip type        |          |                | MTR-85H-PW     |                | MTR-85H-PW     |                |
| - Chart paper                  |          |                | RP-85-PW       |                | RP-85-PW       |                |
| - Ink pen                      |          |                | PG-R-PW        |                | DF-38FP-PW     |                |
| - Recorder housing             |          |                | MDF-S3085-PW   |                | MDF-S3085-PW   |                |
| Small inner door kit           | set of 5 | -              | MDF-5ID5-PW    | MDF-7ID5-PW    | -              |                |
|                                | set of 4 | -              | MDF-5ID4-PW    | MDF-7ID4-PW    | -              |                |

# FROSTLESS FREEZERS

Cost-saving, environmentally friendly sample storage with reduced frost accumulation.

The MDF-DU500ZH-PE and MDF-DU700ZH-PE FrostLess -86°C Upright Freezers with reduced frost build-up on the inner doors. Provides optimum sample storage capacity within a compact footprint combined with natural refrigerants to minimise energy consumption, reduce environmental impact and save money.

## FROSTLESS TECHNOLOGY

The improved packing, insulation materials, and the heated pipe frame of the inner door gasket result in reduced manual labor for freezer maintenance. Less ice will build up on the inner doors and the need for defrosting the freezer is reduced.

## INVERTER COMPRESSORS

While conventional freezers use single speed compressors which cycle on and off, the MDF-DU500ZH and MDF-DU700ZH freezers contain inverter compressors that can run at different speeds to maximise cooling performance under different conditions. Combined with hydrocarbon refrigerants, these compressors ensure the most efficient energy use and reduced heat output.

## INNOVATIVE CABINET DESIGN

The enhanced cabinet design with chamfered edges reduces footprint for use in multi freezer laboratories.



## FROSTLESS TECHNOLOGY



## REDUCING FROST GROWT

The new FrostLess freezers from PHCbi minimise frost build up and with that the time needed for regular maintenance. Reducing frost buildup is critical for smooth running of the freezer, accessibility of the samples and preventing damage of samples.

## NEW THERMAL INSULATION

The new design of the inner doors with VIP PLUS panels and the new thermal insulation technology of the frame suppresses frost growth.



FrostLess Freezers

| Model Number                |          | MDF-DU500ZH-PE   | MDF-DU700ZH-PE    |
|-----------------------------|----------|------------------|-------------------|
| Temperature control range   | °C       |                  | -40 ~ -86         |
| External dimensions (WxDxH) | mm       | 790 x 882 x 1993 | 1030 x 882 x 1993 |
| Internal dimensions (WxDxH) | mm       | 630 x 607 x 1400 | 870 x 607 x 1400  |
| Volume                      | litres   | 525              | 725               |
| Capacity                    | 2" boxes | 352              | 528               |
| Power Consumption           | kWh/day  | 6,1*             | 6,9*              |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load].

# Optional accessories

OPTIONS

FrostLess Freezers

| Model number                   | MDF-DU500ZH-PE | MDF-DU700ZH-PE            |
|--------------------------------|----------------|---------------------------|
| Liquid CO <sub>2</sub> back-up |                | MDF-UB7-PW                |
| Temperature recorders          |                |                           |
| - Circular type                |                | MTR-G85C-PE <sup>1)</sup> |
| - Chart paper                  |                | RP-G85-PW                 |
| - Ink pen                      |                | PG-R-PW                   |
| - Continuous strip type        |                | MTR-85H-PW <sup>1)</sup>  |
| - Chart paper                  |                | RP-85-PW                  |
| - Ink pen                      |                | DF-38FP-PW                |
| - Recorder housing             |                | MDF-S3085-PW              |

1) Requires sensor cover MTR-DU700SF-PW.

# VIP ECO <sup>SMART</sup> FREEZERS

Upgraded with additional features for improved usability.

The VIP ECO SMART -86°C to -40°C upright freezer, along with conventional Password Protection, now offers Face Recognition, and NFC Authentication, coupled with an Electric Lock for enhanced access control. Energy-saving is also further improved by using natural refrigerants combined with inverter compressors.

## Controlled Freezer Access



electronic door lock

Know who entered the freezer and when. Track how long the door was open.

4 Ways to Lock

- Facial recognition
- Keypad
- Proximity card (Optional)
- Unique individual key




facial recognition



## NATURAL REFRIGERANTS

Natural refrigerants, compressors and integrated electronics combined to lower operating costs resulting in less energy consumption and less heat output. Freezer operation is managed by effectively balancing temperature performance and energy management.



## SMART PERFORMANCE

Inverter driven compressors with SMART control allow for superior uniformity by running compressors at slower speeds under steady state conditions. Proprietary algorithms give the control that allows extremely tight temperature uniformity.

## INDUSTRY'S LOWEST ENERGY CONSUMPTION

This freezer delivers lower energy consumption without compromising performance. The VIP ECO SMART uses 30% less energy compared with existing models.

- 5.0 kWh/day, at -80°C
- 4.0 kWh/day, at -70°C



| VIP ECO SMART Freezers      |          |                  |                   |
|-----------------------------|----------|------------------|-------------------|
| Model Number                |          | MDF-DU503VH-PE   | MDF-DU703VH-PE    |
| Temperature control range   | °C       | -40 ~ -86        |                   |
| External dimensions (WxDxH) | mm       | 790 x 882 x 1993 | 1030 x 882 x 1993 |
| Internal dimensions (WxDxH) | mm       | 630 x 600 x 1400 | 870 x 607 x 1400  |
| Volume                      | litres   | 528              | 725               |
| Capacity                    | 2" boxes | 384              | 528               |
| Power Consumption           | kWh/day  | 5,0*             | 5,4*              |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load].

## Optional accessories

| OPTIONS                        |          |                 |                           |
|--------------------------------|----------|-----------------|---------------------------|
| VIP ECO SMART Freezers         |          |                 |                           |
| Model number                   |          | MDF-DU503VH-PEE | MDF-DU703VH-PE            |
| Small inner door kit           | set of 5 |                 | MDF-5ID5-PW <sup>1)</sup> |
|                                | set of 4 |                 | MDF-5ID4-PW               |
| NFC System <sup>2)</sup>       |          |                 | NFC-reader: MTR-NFC-PW    |
|                                | set of 5 |                 | NFC-Cards: MTR-IC-PW      |
| Liquid CO <sub>2</sub> back-up |          |                 | MDF-UB9-PW                |
| Temperature recorders          |          |                 |                           |
| - Circular type                |          |                 | MTR-G85C-PE <sup>3)</sup> |
| - Chart paper                  |          |                 | RP-G85-PW                 |
| - Ink pen                      |          |                 | PG-R-PW                   |
| - Continuous strip type        |          |                 | MTR-85H-PW <sup>3)</sup>  |
| - Chart paper                  |          |                 | RP-85-PW                  |
| - Ink pen                      |          |                 | DF-38FP-PW                |
| - Recorder housing             |          |                 | MDF-S3085-PW              |
| Recorder sensor cover          |          |                 | MTR-DU700SF-PW            |

1) Installation of small inner door kit may effect usable storage capacity.

2) The NFC System is only available in EEA member countries, UK, Switzerland and Turkey.

3) Requires sensor cover MTR-DU700SF-PW. handle and other external projections.

# VIP ECO ULT FREEZERS

Cost-saving and environmentally friendly sample storage within an optimal footprint.

VIP ECO Ultra Low Temperature Freezers with natural refrigerants minimise energy consumption, reduce environmental impact and save money. Innovative technology provide secure storage of valuable research and clinical samples. The VIP PLUS vacuum insulation ensures an optimal footprint to storage capacity ratio.

Leveraging the power of natural hydrocarbon refrigerants also allows the VIP ECO ULT Freezers to use smaller compressors, and reduce energy consumption. The natural hydrocarbon refrigerants combined with VIP insulation technology also help the environment by reducing the carbon footprint with up to 40% fewer emissions.

## Inverter Compressors



The intelligent control of the inverter compressor optimizes running speed. When the inverter compressor is running as normal it will stay on for longer than a conventional compressor but at a minimal speed. This reduces the power consumption and keeps freezer temperatures stable.

## Natural Refrigerants



Naturally occurring hydrocarbon (HC) refrigerants improve performance and reduce running costs.

When low electrical running costs and environmental considerations are of paramount importance, put your trust in VIP ECO and PRO ECO ULT freezers.

- New heat exchanger design for greater surface area contact and overall efficiency.
- Advanced space-saving VIP PLUS Insulation (compared to our conventional models).
- Graphical LCD display with data monitoring and data log exported by USB (for MDF-DU502VH-PE, MDF-DU502VHW-PE, MDF-DU702VH-PE and MDF-DU702VHW-PE).
- VIP Plus Insulation
- Low heat emission
- Two Vacuum release ports (MDF-DU901VHL-PE: 1 Vacuum release port)
- No icing on frame
- Alarm and Safety functions





## VIP PLUS INSULATION



PHCbi's patented VIP PLUS technology has resulted in a revolutionary vacuum insulation cabinet construction with improved thermal properties for superior temperature performance.

The VIP ECO Series also uses vacuum insulation panel (VIP) technology reducing wall thickness by around 50%, achieving 30% more storage capacity, and reducing the average cost per box stored. Leveraging the power of natural hydrocarbon refrigerants also allows the VIP ECO ULT freezers to use smaller compressors, due to their greater efficiency. The natural hydrocarbon refrigerants combined with VIP PLUS insulation technology also help the environment by reducing the carbon footprint with up to 40% fewer emissions.

## REDUCED RUNNING COSTS

The use of highly efficient hydrocarbon refrigerants results in reduced energy consumption and lower running costs. With key equipment and instrumentation operating continuously laboratories are able to significantly reduce running costs by investing in energy efficient facilities. PHC designs and builds advanced preservation systems to deliver maximum cost efficiency while maintaining the reliability and performance necessary for reliable storage of valuable research and clinical samples.


## EXTREMELY LOW ENVIRONMENTAL IMPACT

Naturally occurring hydrocarbon (HC) refrigerants used within the VIP ECO ULT Freezers are non ozone depleting, have short atmospheric lifetimes and have extremely low global warming potentials (GWP's). This makes the freezers very environmentally friendly so they are an ideal solution for complying with objectives for reduced carbon footprints.

## OPTIMUM UNIFORMITY

Uneven interior temperatures can lead to a loss in sample integrity. PHCbi freezers with uniform, stable temperatures and quick recovery times provide the best protection for your samples, ensuring reliable preservation while guarding against degradation.

Surpasses the customer preference of **+/-5°C\***




\* Based on internal validation data tested at -80°C setpoint, in an empty chamber with 23°C ambient temperature.

\* The data may vary depending on the use, circumstances and optional accessories. Validation documents can be provided for each serial number for an additional fee.

MDF-DU502VH up to **55%** more efficient  
MDF-DU702VH up to **46%** more efficient

- 1) Model MDF-DU702VH-PE (728 litres) compared to MDF-U73V at set value -80°C, 23°C ambient, no load, 230V 50Hz, € 0.12/kWh. Actual energy use and savings will depend on operating conditions and price of electricity paid.
- 2) Model MDF-DU502VH-PE (526 litres) compared to MDF-U53V at set value -80°C, 23°C ambient, no load, 230V 50Hz, € 0.12/kWh. Actual energy use and savings will depend on operating conditions and price of electricity paid.



## VACUUM RELEASE PORT

A vacuum release port allows smooth door opening when the door seal is tightened by negative pressure generated from temperature difference between chamber and ambient.



## FLEXIBLE SHELF LAYOUT

Multiple shelf configurations in the upright models allow a variety of storage options. Organize your samples by using your existing inventory racks or select from the many different rack types we offer. PHCbi's racks are made of stainless steel or anodized aluminum. The aluminum racks are very light, yet sturdy and corrosion free.

## EZlatch



The EZ Latch door handle was developed based on human engineering. It is designed for secure door opening/closing with Minimum effort as the name suggests and make access to stored sample easier.

## MEETING YOUR FREEZER STORAGE NEEDS

### An organized freezer will provide you with:

- Time savings - locate, retrieve and replace your samples easily and quickly.
- Cost savings - organized samples and cell lines can help to reduce the number of freezers.
- Added sample security and energy savings - samples are better protected and are less exposed to ambient temperatures as door opening times can be reduced when placing and retrieving samples, which also reduces energy use.

For an overview of the racks designed for the VIP ECO series see page 42-46.

PHCbi racks\* are made of stainless steel or anodized aluminum. The aluminum racks are very light, yet sturdy and corrosion free.


### Meeting your freezer storage needs

- High quality racks – designed for safe working and easy access to samples.
- Affordable solutions – making freezer storage cost-effective as well as space-efficient.
- Large selection of products – additional rack types and boxes are available on request.



HCS-296

NIR-224U





VIP ECO ULT Freezers

| Model Number                |          | MDF-DU502VH-PE   | MDF-DU702VH-PE    | MDF-DU901VHL-PE   |
|-----------------------------|----------|------------------|-------------------|-------------------|
| Temperature control range   | °C       | -40 ~-86         |                   | -50 ~-86          |
| External dimensions (WxDxH) | mm       | 790 x 882 x 1993 | 1030 x 882 x 1993 | 1150 x 870 x 1993 |
| Internal dimensions (WxDxH) | mm       | 630 x 600 x 1400 | 870 x 600 x 1400  | 1010 x 600 x 1400 |
| Volume                      | litres   | 528              | 729               | 845               |
| Capacity                    | 2" boxes | 384              | 576               | 672               |
| Power Consumption           | kWh/day  | 6,7*             | 7,7*              | 8,7*              |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load].

## Optional accessories

OPTIONS

VIP ECO ULT Freezers

| Model number                   |          | MDF-DU502VH-PE | MDF-DU702VH-PE          | MDF-DU901VHL-PE    |
|--------------------------------|----------|----------------|-------------------------|--------------------|
| Liquid CO <sub>2</sub> back-up |          |                | MDF-UB7-PW              |                    |
| Temperature recorders          |          |                |                         |                    |
| - Circular type                |          |                | MTR-G85C-PE             |                    |
| - Chart paper                  |          |                | RP-G85-PW <sup>8)</sup> |                    |
| - Ink pen                      |          |                | PG-R-PW                 |                    |
| - Continuous strip type        |          |                | MTR-85H-PW              |                    |
| - Chart paper                  |          |                | RP-85-PW <sup>8)</sup>  |                    |
| - Ink pen                      |          |                | DF-38FP-PW              |                    |
| - Recorder housing             |          |                | MDF-S3085-PW            |                    |
| Drawers                        | qty      | -              | -                       | -                  |
| Small inner door kit           | set of 2 | -              | -                       | MDF-9ID-PW (max 2) |
|                                | set of 5 | MDF-5ID5-PW    | MDF-7ID5-PW             | -                  |
|                                | set of 4 | MDF-5ID4-PW    | MDF-7ID4-PW             | -                  |

# VIP ECO / PRO ECO ULT FREEZERS

These ultra low temperature freezers offer advanced cabinet design, reliable refrigeration systems and easy-to-use controllers making them ideal for the long-term secure storage. Every component is carefully selected and matched for optimum operation under demanding laboratory conditions, while the internal layout of the refrigeration system is meticulously designed for maximum heat removal, reducing stress on the system and therefore providing the highest levels of reliability and durability.

When sample security and peace of mind are of paramount importance, put your trust in VIP PLUS ultra low freezers.

- Personal sized ult chest freezer
- User-friendly filterless design
- Quiet operation
- Excellent sample security
- Uniform sample storage
- VIP PLUS insulation maximizes storage capacity

Add minimising energy consumption and reducing environmental impact to sample security, you get PRO ECO ULT Freezers.

- Energy efficiency
- Lower running costs
- Reduced environmental impact
- Excellent sample security
- Natural Refrigerants
- Heat Exchanger Design

## VIP PLUS INSULATION



PHCbi's patented VIP PLUS technology has resulted in a revolutionary vacuum insulation cabinet construction with improved thermal properties for superior temperature performance.



Model: MDF-DC102VH-PE

## Natural Refrigerants

Naturally occurring hydrocarbon (HC) refrigerants improve performance and reduce running costs.



Model: MDF-DU300H-PE



|                             |          | VIP ECO Freezers |                  | PRO ECO ULT Freezers |
|-----------------------------|----------|------------------|------------------|----------------------|
| Model Number                |          | MDF-DC102VH-PE   | MDF-DC202VJ-PE   | MDF-DU300H-PE        |
| Temperature control range   | °C       | 558 x 688 x 993  | -40 to -86       | -86                  |
| External dimensions (WxDxH) | mm       | 405 x 490 x 426  | 1018 x 688 x 993 | 750 x 870 x 1830     |
| Internal dimensions (WxDxH) | mm       | 84               | 865 x 490 x 426  | 490 x 600 x 1140     |
| Volume                      | litres   | 42               | 180              | 333                  |
| Capacity                    | 2" boxes | 4,6*             | 105              | 216                  |
| Power Consumption           | kWh/day  |                  | 6,4*             | 5,7*                 |

\* [Set value temp. -80°C, Ambient temp. 23 °C, no load]

## MDF-DC102VH-PE & MDF-DC202VH-PE

PHCbi's MDF-DC102VH-PE and MDF-DC202VH-PE VIP ECO Chest Freezers provide maximum sample storage capacity within an optimum footprint combined with industry leading cooling performance and reliability for an ideal long-term preservation solution.

## MDF-DU300H-PE

The PRO ECO -86°C Ultra-low Freezer, MDF-DU300H-PE, is part of our range ECO ultra-low freezers with natural refrigerants. They minimise energy consumption, reduce environmental impact and save money. Innovative technology and Class IIa medical device certification provide secure storage of valuable research and clinical samples.

# Optional accessories

|                                |     | VIP ECO Freezers |                | PRO ECO ULT Freezers |
|--------------------------------|-----|------------------|----------------|----------------------|
| Model Number                   |     | MDF-DC102VH-PE   | MDF-DC202VH-PE | MDF-DU300H-PE        |
| Liquid CO <sub>2</sub> back-up |     |                  | MDF-UB8-PE     | CVK-UB2-PW           |
| Temperature recorders          |     |                  |                |                      |
| - Circular type <sup>1)</sup>  |     |                  | MTR-G85C-PE    | MTR-G85C-PE          |
| - Chart paper                  |     |                  | RP-G85-PW      | RP-G85-PW            |
| - Ink pen                      |     |                  | PG-R-PW        | PG-R-PW              |
| - Continuous strip type        |     |                  | MTR-85H-PW     | MTR-85H-PW           |
| - Chart paper                  |     |                  | RP-85-PW       | RP-85-PW             |
| - Ink pen                      |     |                  | DF-38FP-PW     | DF-38FP-PW           |
| - Recorder housing             |     |                  | MDF-S3085-PW   | MDF-S3085-PW         |
| Drawers                        | qty |                  | -              | MDF30RPW (max)2      |

<sup>1)</sup> Requires sensor cover MTR-C8-PW.

# HYBRID ULT FREEZERS

When a high-quality ultra low temperature or cryogenic freezer is equipped with a Hybrid water cooling option, the unit can handle a chilled water circuit to extract the generated heat from the condenser or use the traditional air-cooled way with a fan motor. With this new setup, the freezer can switch from water cooled to air cooled in case the water system is not operated. A Hybrid water cooled freezer will contribute to a significant reduction in power consumption and will also reduce the amount of heat dissipated into the air. Compared to an air cooled freezer a Hybrid water cooled freezer will also have an improved temperature stability.

When reduced running costs and environmental objectives are of paramount importance, put your trust in HYBRID ULT freezers.


- Lower energy use
- Cost saving
- Re-use of energy
- Reduced heat dissipation
- Improved sample protection
- Faster recovery time
- Reduction of airconditioning

## HYBRID TECHNOLOGY



PHCbi's Hybrid water cooled technology on ECO VIP ultra low temperature and cryogenic freezers improves the compressor efficiency.

The power consumption of a Hybrid ULT can be reduced by typically up to 14% compared to the equivalent air-cooled model. PHC offers a selection of VIP ECO and cryogenic freezers with hybrid water cooling technology.




Model: MDF-DU702VHW-PE



Model: MDF-DU502VHW-PE

## LOWER ENERGY USE AND RUNNING COSTS

- Water has a greater heat absorption capacity than air. Therefore a water cooled condenser is more efficient than a conventional air cooled condenser. As a result, the compressor efficiency is improved, and the power consumption of a Hybrid ULT can be reduced by typically up to 14% compared to the equivalent air cooled model.
- Compared to an air cooled freezer a Hybrid water cooled freezer will reduce the heat dissipation into the air. This results in the reduction of air conditioning requirements for further cost savings.
- Potential to re-use the energy in the Hybrid water cooling system for other purposes within the facility to reduce site-wide energy costs.



## IMPROVED FREEZER PERFORMANCE AND SAMPLE PROTECTION


- The greater cooling capacity of water compared to air improves the performance of the refrigeration system leading to reduced pull-down times. This provides



Model: MDF-C2156VANW-PE

faster temperature recovery after door opening and sample loading, thereby protecting samples by helping to maintain the correct temperature.


- If the room air-conditioning fails, the room with Hybrid water cooled freezers will not warm up so quickly so samples in the freezers will remain safe for longer.
- The installation of a water cooled system to remove heat from ultra-low temperature and cryogenic freezers can assist organisations to reduce running costs and to meet environmental and energy reduction objectives.
- If the water system stops for any reason (failure or maintenance) the Hybrid water cooled freezer will start to work automatically on the air cooled condenser. This provides the best protection for your samples.
- The inverter compressors inside the MDF-DU502VHW-PE and MDF-DU702VHW-PE are developed for the best performance with the lowest energy consumption.



## PLATE HEAT EXCHANGER

Heat energy from inside the freezer compartment is transferred by refrigerant gasses to a plate heat exchanger. Inside the plate heat exchanger, energy is transferred from the refrigerant to a closed water circuit. The greater cooling capacity of water compared to air improves the performance of the refrigeration system leading to reduced pull-down times. This provides faster temperature recovery after door opening and sample loading.

## HYBRID CASCADE COOLING SYSTEM



## REDUCED HEAT DISSIPATION

### COMPARISON MDF-DU702VH-PE VS MDF-DU702VHW-PE

| Model  | Heat dissipation  | Power consumption (kWh/day) | Power consumption (W/hour) | Heat dissipation (kcal/hour) | Percentage (%) |
|--|-------------------|-----------------------------|----------------------------|------------------------------|----------------|
| MDF-DU702VH-PE   | Total (air)       | 9.4                         | 391.7                      | 336.8                        | 100%           |
| MDF-DU702VHW-PE  | Total (air+water) | 7.8                         | 325.0                      | 279.4                        | 100%           |
|  | Air               |                             |                            | 105.4                        | 38%            |
|  | Water             |                             |                            | 174.0                        | 62%            |
| Reduction of Heat Dissipation into the air (MDF-DU702VH-PE vs MDF-DU702VHW-PE) |                   |                             |                            |                              | 69%            |
| Reduction of Power Consumption (MDF-DU702VH-PE vs MDF-DU702VHW-PE)             |                   |                             |                            |                              | 17%            |

Data measured with freezers running at Setpoint -80°C, Ambient Temperature 27.2°C, No Load inside the freezers

Conclusion: The hybrid cooling system of a MDF-DU702VHW-PE contributes in the reduction of Power Consumption and Heat Dissipation to the air compared to a standard air cooled MDF-DU702VH-PE. Especially at high ambient temperatures like 27.2 °C the reduction in Power Consumption and Heat Dissipation into the air is significant.





| VIP HYBRID ULT Freezers     |          |                  |                   |
|-----------------------------|----------|------------------|-------------------|
| Model Number                |          | MDF-DU502VHW-PE  | MDF-DU702VHW-PE   |
| Temperature control range   | °C       | -40 ~ -86        |                   |
| External dimensions (WxDxH) | mm       | 790 x 882 x 1993 | 1030 x 882 x 1993 |
| Internal dimensions (WxDxH) | mm       | 630 x 600 x 1400 | 870 x 600 x 1400  |
| Volume                      | litres   | 528              | 729               |
| Capacity                    | 2" boxes | 384              | 576               |

| Cryogenic Hybrid Freezer    |          |                   |
|-----------------------------|----------|-------------------|
| Model Number                |          | MDF-C2156VANW-PE  |
| Temperature control range   | °C       | -150              |
| External dimensions (WxDxH) | mm       | 1730 x 765 x 1010 |
| Internal dimensions (WxDxH) | mm       | 760 x 495 x 615   |
| Volume                      | litres   | 231               |
| Capacity                    | 2" boxes | 150               |
| Net weight (approx)         | kg       | 318               |



## Optional accessories

### OPTIONS

| Hybrid ULT Freezers            |          |                 |                         |                      |
|--------------------------------|----------|-----------------|-------------------------|----------------------|
| Model number                   |          | MDF-DU502VHW-PE | MDF-DU702VHW-PE         | MDF-C2156VANW-PE     |
| Liquid CO <sub>2</sub> back-up |          |                 | MDF-UB7-PW              | -                    |
| Liquid N <sub>2</sub> back-up  |          |                 | -                       | Supplied as standard |
| Temperature recorders          |          |                 |                         |                      |
| - Circular type                |          |                 | MTR-G85C-PE             | MTR-155H-PW          |
| - Chart paper                  |          |                 | RP-G85-PW <sup>8)</sup> | RP-155-PW            |
| - Ink pen                      |          |                 | PG-R-PW                 | DF-38FP-PW           |
| - Continuous strip type        |          |                 | MTR-85H-PW              | MDF-S30150-PW        |
| - Chart paper                  |          |                 | RP-85-PW <sup>8)</sup>  |                      |
| - Ink pen                      |          |                 | DF-38FP-PW              |                      |
| - Recorder housing             |          |                 | MDF-S3085-PW            |                      |
| Small inner door kit           | set of 5 | MDF-5ID5-PW     | MDF-7ID5-PW             |                      |
|                                | set of 4 | MDF-5ID4-PW     | MDF-7ID4-PW             |                      |

# CRYOGENIC VIALS

Cryogenic Vials for cryopreservation and transport can be used to safely store valuable samples such as cells, serum, and antibodies.

The cryogenic vials have a silicone packing which is attached to the cap of the vials to ensure high airtightness and prevent drying of samples and liquid leakage.

## FEATURES OF CRYOGENIC VIALS

- Excellent airtightness, leak tested in accordance with IATA (International Air Transport Association) international standards, for air transportation packaging.
- Can be stored for a prolonged period of time at -190 °C.
- Easy loading and unloading of samples prevent contamination of the outer wall of the tube due to dripping.
- Various special racks for easy handling are also available.

## AVAILABLE PRODUCTS

- 2 types of caps: outer cap and inner cap
- 4 capacities inner cap: 1.2 ml, 2 ml, 4 ml, 5 ml
- 4 capacities outer cap: 1 ml, 2 ml, 4 ml, 5 ml
- 5 cap colours: white, red, yellow, blue, green (a 5-colour mix set is also available)

### Various special boxes are also available

- MS-7550L, MS-7550G, MS-7550P, MS-75600

### Lock Stand: MS-78000

- Locked at its bottom, the tube's cap can be opened and closed with one hand.

## SLIM TUBES

Our original "Slim Tube" is a **centrifuge compatible product**.

Using regular tubes you need to transfer samples to a different tube which is centrifuge compatible, this increases the risk of contamination whilst transferring the samples. Since the slim tube is a centrifuge compatible tube, the risk of contamination due to transferring tubes is avoided. Centrifugal operation is possible due to the narrow body diameter. Max Centrifuge force: 19,500 G



## PRODUCT TYPE



- Outer cap type



- Inner cap type

## PRODUCT SPECIFICATIONS

- Material: polypropylene.
- Sterilization: radiation sterilized.
- Pyrogen & endotoxin Free: endotoxin level less than 0.5 EU/mL (in accordance with the Japanese Pharmacopoeia endotoxin test method).
- Sealing pressure: 95kPa pressure tested (In accordance with IATA PI602/650). Leak test conducted by a third-party organization.
- Liquid nitrogen: storage in gas phase possible.
- Autoclave: yes (121°C, 20 min).
- Cytotoxicity: no toxicity (Intraocular lens approval standard. Cell growth inhibition test. Agar diffusion test).
- Mutagenicity: no inducibility (Standard for mutagenicity test using microorganisms).
- Exothermic substances: compliant (Exothermic substances test).

| Outer cap type | Product code | Product name        | Shape         |                               | Quantity              |
|----------------|--------------|---------------------|---------------|-------------------------------|-----------------------|
|                |              |                     | Bottom        | Size                          |                       |
|                | MS-4601      | Cryogenic vial 1 ml | Self-standing | Outer diameter Ø 12.5 x 43 mm | 500 (50 pcs/bag x 10) |
|                | MS-4603      | Cryogenic vial 2 ml | Self-standing | Outer diameter Ø 12.5 x 46 mm | 500 (50 pcs/bag x 10) |
|                | MS-4604      | Cryogenic vial 4 ml | Self-standing | Outer diameter Ø 12.5 x 71 mm | 300 (30 pcs/bag x 10) |
|                | MS-4605      | Cryogenic vial 5 ml | Self-standing | Outer diameter Ø 12.5 x 87 mm | 300 (30 pcs/bag x 10) |

Cap colour code (W: white, R: red, Y: yellow, B: blue, G: green, X:5 colour set). In ordering, please indicate colours after product code (e.g. MS-4601W, MS-4601R, etc.).

| Inner cap type | Product code | Product name          | Shape         |                               | Quantity              |
|----------------|--------------|-----------------------|---------------|-------------------------------|-----------------------|
|                |              |                       | Bottom        | Size                          |                       |
|                | MS-4501      | Cryogenic vial 1.2 ml | Self-standing | Outer diameter Ø 12.5 x 42 mm | 500 (50 pcs/bag x 10) |
|                | MS-4502      | Cryogenic vial 2 ml   | Round         | Outer diameter Ø 12.5 x 48 mm | 500 (50 pcs/bag x 10) |
|                | MS-4503      | Cryogenic vial 2 ml   | Self-standing | Outer diameter Ø 12.5 x 48 mm | 500 (50 pcs/bag x 10) |
|                | MS-4504      | Cryogenic vial 4 ml   | Round         | Outer diameter Ø 12.5 x 70 mm | 300 (30 pcs/bag x 10) |
|                | MS-4505      | Cryogenic vial 5 ml   | Round         | Outer diameter Ø 12.5 x 90 mm | 300 (30 pcs/bag x 10) |

Cap colour code (W: white, R: red, Y: yellow, B: blue, G: green, X:5 colour set). In ordering, please indicate colours after product code (e.g. MS-4501W, MS-4501R, etc.).

| Slim tube type | Product code | Product name     | Shape         |                               | Quantity              |
|----------------|--------------|------------------|---------------|-------------------------------|-----------------------|
|                |              |                  | Bottom        | Size                          |                       |
|                | MS-4701      | Slim Tube 0.5 ml | Self-standing | Outer diameter Ø 10.5 x 46 mm | 500 (50 pcs/bag x 10) |
|                | MS-4702      | Slim Tube 1.5 ml | Self-standing | Outer diameter Ø 10.5 x 46 mm | 500 (50 pcs/bag x 10) |

| Cryogenic Vials Accessories | Product code                      | Product name  | Product description  | Quantity            |
|-----------------------------|-----------------------------------|---|--|---------------------|
|                             | MS-7550LZ                         | Rack for Cryogenic vials with lid   | Colour: grey   | 10 (10 pcs/box x 1) |
|                             | MS-7550GZ                         | Rack for Cryogenic vials with lid   | Colour: green  | 10 (10 pcs/box x 1) |
|                             | MS-7550PZ                         | Rack for Cryogenic vials with lid   | Colour: pink   | 10 (10 pcs/box x 1) |
|                             | MS-75600Z                         | Rack for Cryogenic vials with lid   | Colour: white  | 10 (10 pcs/box x 1) |
|                             | MS-99150Z                         | Box for Cryogenic vials LH  | <ul style="list-style-type: none"> <li>LH: L (Large size), H (High) • 146(W) x 146(L) x 76(H) mm,</li> <li>for 4mL &amp; 5mL vials with capacity of 100 vials</li> </ul> | 20 (20 pcs/box x 1) |
|                             | MS-99160Z                         | Box for Cryogenic vials LL  | <ul style="list-style-type: none"> <li>LL: L (Large size), L (Low) • 146(w) x 146(L) x 51(H) mm</li> <li>for 1mL &amp; 2mL vials with capacity of 100 vials</li> </ul>   | 20 (20 pcs/box x 1) |
|                             | MS-99170Z                         | Box for Cryogenic vials M   | <ul style="list-style-type: none"> <li>M: Medium • 133(w) x 133(L) x 51(H) mm</li> <li>for 1mL &amp; 2mL vials with capacity of 81 vials</li> </ul>                      | 20 (20 pcs/box x 1) |
| MS-75200Z                   | Coloured caps for Cryogenic vials | <ul style="list-style-type: none"> <li>8 colour caps for organized storage of inner cap typed vials</li> <li>Snap-in cap head for easy identification</li> <li>Making it possible to be identified on the surface due to the different colours</li> <li>Cannot be autoclaved</li> </ul> | 1600 (1600 pcs/box x 1)  |                     |
| MS-78000Z                   | Lock Stand for Cryogenic vials    | <ul style="list-style-type: none"> <li>Has a locking mechanism which allows one-handed operation to open and close the cap of vials</li> <li>Autoclavable (121°C, 20 minutes)</li> </ul>  | 5 (5 pcs/box x 1)  |                     |

# INVENTORY RACKS

Optimize and simplify access to stored materials with inventory racks suitable for your application.

Organizing the space in your freezer can help you to make your work more efficient. Good sample management can result in both cost and time savings. By choosing the right racks, not only will efficiency be increased, but the risk of sample degradation due to exposure to ambient temperatures will be significantly reduced. This is a fact for any type of cell or tissue sample.

## DURABLE LONG-TERM STORAGE

Racks are available in stainless steel and aluminium. Both racks types provide a durable and reliable long-term storage solution.

## EASY SAMPLE ACCESS

Intuitive rack and handle design



## MAXIMIZE FREEZER BOX STORAGE

Optimum freezer and footprint cost per box

## EXCLUSIVELY DESIGNED FOR PHCbi

Made for upright and chest freezers

- Stainless steel or aluminium construction designed to withstand a lifetime of use.
- Racks are available for use with fibreboard or plastic boxes.
- Intuitive rack designs make it easy to access and return racks to the freezer.
- Maximum storage in freezer footprint.
- Sized to fit the PHCbi family of upright and chest freezers.
- Custom made racks to suit your unique storage requirements are available upon request.

### Boxes and dividers

- Standard 2" and 3" boxes are available, in sturdy moisture resistant cardboard and in polypropylene boxes
- For the XC series cryosystems we have the smaller 2" minibox including a 25 cell divider (B2CM + D25M).
- Sturdy cardboard cell dividers come in a wide selection of sizes to accommodate the storage of a variety of tubes and vials.



## CRYOGENIC INVENTORY RACKS

Models: CBS

### STANDARD SQUARE RACKS

| Model number  | Rack type  | Description  | Rack/quantity                                       |
|---------------|------------|--|---|
| V-1500AB-ST   | 2001A-100S | Aluminium rack system with cardboard boxes + dividers. | 7 racks x 13 boxes high. Max. cap. 9.100 2ml vials. |
| S-1500AB-ST   | 2001S-C81  | Aluminium rack system with cardboard boxes + dividers. | 7 racks x 13 boxes high. Max. cap. 9.100 2ml vials. |
| V-3000AB-ST   | 3101A-100S | Aluminium rack system with cardboard boxes + dividers. | 17 racks x 13 boxes high. cap. 22.100 2ml vials.    |
| S-3000AB-ST   | 3101A-100S | Aluminium rack system with cardboard boxes + dividers. | 17 racks x 13 boxes high. cap. 22.100 2ml vials.    |
| V-5000AB-ST   | 3301A-100S | Aluminium rack system with cardboard boxes + dividers. | 28 racks x 13 boxes high. cap. 36.400 2ml vials.    |
| S-5000AB-ST   | 3301A-100S | Aluminium rack system with cardboard boxes + dividers. | 28 racks x 13 boxes high. cap. 36.400 2ml vials.    |
| V-5000ABEH-ST | 3325A-100S | Aluminium rack system with cardboard boxes + dividers. | 28 racks x 15 boxes high. cap. 42.000 2ml vials.    |
| S-5000ABEH-ST | 3325A-100S | Aluminium rack system with cardboard boxes + dividers. | 28 racks x 15 boxes high. cap. 42.000 2ml vials.    |

### VERTICAL RACKS

| Model number  | Rack type           | Description  | Rack/quantity                                   |
|---------------|---------------------|--|---|
| V-1500AB-ST   | RC-V1500-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 20 racks x 5 boxes high. cap. 10.000 2ml vials. |
| S-1500AB-ST   | RC-S1500-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 20 racks x 5 boxes high. cap. 10.000 2ml vials  |
| V-3000AB-ST   | RC-V3000-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 48 racks x 5 boxes high. cap. 24.000 2ml vials. |
| S-3000AB-ST   | RC-S3000-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 48 racks x 5 boxes high. cap. 24.000 2ml vials. |
| V-5000AB-ST   | RC-V5000-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 80 racks x 5 boxes high. cap. 40.000 2ml vials. |
| S-5000AB-ST   | RC-S5000-1209-VLR   | Stainless steel rack system with cardboard boxes + dividers. | 80 racks x 5 boxes high. cap. 40.000 2ml vials  |
| V-5000ABEH-ST | RC-V5000EH-1208-VLR | Stainless steel rack system with cardboard boxes + dividers. | 80 racks x 6 boxes high. cap. 48.000 2ml vials. |
| S-5000ABEH-ST | RC-S5000EH-1208-VLR | Stainless steel rack system with cardboard boxes + dividers. | 80 racks x 6 boxes high. cap. 48.000 2ml vials  |

## CRYOGENIC MDF FREEZERS

Models: MDF-1156(ATN)-PE | MDF-C2156VAN-PE

### ALUMINIUM RACK SOLUTIONS

| Vertical rack type      | Box type | Rack/quantity<br>Aluminium | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |        | Maximum box dimensions            |                          |                           |
|-------------------------|----------|----------------------------|------------------------------------|----------------|------------------|-------|--------|-----------------------------------|--------------------------|---------------------------|
|                         |          |                            |                                    |                | Width            | Depth | Height | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| <b>MDF-1156(ATN)-PE</b> |          |                            |                                    |                |                  |       |        |                                   |                          |                           |
| side opening            | 2 inch   | 9 x NIR-209C               | 3x3                                | 81             | 142              | 141   | 516    | 133                               | 133                      | 53                        |
| side opening            | 3 inch   | 9 x NIR-306C               | 3x3                                | 54             | 142              | 141   | 516    | 133                               | 133                      | 75                        |
| <b>MDF-C2156VAN-PE</b>  |          |                            |                                    |                |                  |       |        |                                   |                          |                           |
| opening                 | 2 inch   | 15 x NIR-210C              | 5x3                                | 150            | 142              | 141   | 590    | 133                               | 133                      | 53                        |
| side opening            | 3 inch   | 15 x NIR-307C              | 5x3                                | 105            | 142              | 141   | 590    | 133                               | 133                      | 78                        |

### STAINLESS STEEL RACK SOLUTIONS

| Vertical rack type      | Box type | Rack/quantity<br>Stainless steel | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |        | Maximum box dimensions            |                          |                           |
|-------------------------|----------|----------------------------------|------------------------------------|----------------|------------------|-------|--------|-----------------------------------|--------------------------|---------------------------|
|                         |          |                                  |                                    |                | Width            | Depth | Height | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| <b>MDF-1156(ATN)-PE</b> |          |                                  |                                    |                |                  |       |        |                                   |                          |                           |
| side opening            | 2 inch   | 9 x SCR-102-N                    | 3x3                                | 90             | 139.7            | 144   | 564.13 | 136                               | 142                      | 54                        |
| side opening            | 3 inch   | 9 x SCR-063-N                    | 3x3                                | 54             | 139.7            | 144   | 494.28 | 136                               | 142                      |                           |
| <b>MDF-C2156VAN-PE</b>  |          |                                  |                                    |                |                  |       |        |                                   |                          |                           |
| side opening            | 2 inch   | 15 x SCR-102-N                   | 5x3                                | 150            | 139.7            | 144   | 564.13 | 136                               | 142                      | 54                        |
| side opening            | 3 inch   | 15 x SCR-073-N                   | 5x3                                | 105            | 139.7            | 144   | 575.31 | 136                               | 142                      | 75                        |

\* Unit: mm

# TWINGUARD / VIP ECO INVENTORY RACKS

Models: MDF-DU302VX-PE | MDF-DU502VX-PE | MDF-DU702VX-PE | MDF-DU502VH-PE | MDF-DU702VH-PE | MDF-DU503VH-PE  
MDF-DU703VH-PE | MDF-DU901VHL-PE

## ALUMINIUM RACK SOLUTIONS

| Vertical rack type   | Box type | Rack/quantity<br>Aluminium    | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |         | Maximum box dimensions            |                          |                           |
|--|----------|-------------------------------|------------------------------------|----------------|------------------|-------|---------|-----------------------------------|--------------------------|---------------------------|
|  |          |                               |                                    |                | Width            | Depth | Height  | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| <b>MDF-DU302VX-PE</b>  |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays   | 2 inch   | 12 x HCS-32-5584/143          | 3x4                                | 240            | 143              | 560   | 280     | 133                               | 136                      | 52                        |
| side opening   | 2 inch   | 12 x NIR-220U                 | 3x4                                | 240            | 139              | 559   | 279     | 135                               | 135                      | 52                        |
| with trays   | 3 inch   | 12 x HCS-32-3804/143          | 3x4                                | 144            | 143              | 560   | 232     | 133                               | 136                      | 71                        |
| side opening   | 3 inch   | 12 x NIR-312U                 | 3x4                                | 144            | 139              | 559   | 279     | 135                               | 135                      | 88                        |
| <b>MDF-DU502VX-PE / MDF-DU502VH-PE / MDF-DU502VH-PE / MDF-DU503VH-PE</b> |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays   | 2 inch   | 4 x HCS-296                   | 2x2                                | 384            | 280              | 560   | 685     | 133                               | 136                      | 52                        |
| with trays   | 2 inch   | 16 x HCS-6564                 | 4x4                                | 384            | 140              | 560   | 339     | 130                               | 133                      | 52                        |
| side opening   | 2 inch   | 16 x NIR-224U                 | 4x4                                | 384            | 139              | 559   | 334     | 135                               | 135                      | 52                        |
| with trays   | 3 inch   | 16 x HCS-4804                 | 4x4                                | 256            | 140              | 560   | 320     | 130                               | 133                      | 75                        |
| side opening   | 3 inch   | 16 x NIR-316U                 | 4x4                                | 256            | 139              | 560   | 334     | 135                               | 135                      | 75                        |
| <b>MDF-DU702VX-PE / MDF-DU702VH-PE / MDF-DU702VH-PE / MDF-DU703VH-PE</b> |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays   | 2 inch   | 6 x HCS-296                   | 3x2                                | 576            | 280              | 560   | 685     | 133                               | 136                      | 52                        |
| with trays   | 2 inch   | 24 x HCS-6564                 | 6x4                                | 576            | 140              | 560   | 339     | 130                               | 133                      | 52                        |
| side opening   | 2 inch   | 24 x NIR-224U                 | 6x4                                | 576            | 139              | 559   | 334     | 135                               | 135                      | 52                        |
| with trays   | 3 inch   | 24 x HCS-4804                 | 6x4                                | 384            | 140              | 560   | 320     | 130                               | 133                      | 75                        |
| side opening   | 3 inch   | 24 x NIR-316U                 | 6x4                                | 384            | 139              | 559   | 324     | 135                               | 135                      | 75                        |
| <b>MDF-U901VHL-PE</b>  |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays   | 2 inch   | 14 x HCS-5584 + 14 x HCS-6564 | 7x2 + 7x2                          | 616            | 140              | 560   | 290/339 | 130                               | 133                      | 52                        |
| side opening   | 2 inch   | 14 x NIR-220U + 14 x NIR-224U | 7x2 + 7x2                          | 616            | 139              | 559   | 279/334 | 130                               | 133                      | 52                        |
| with trays   | 3 inch   | 14 x HCS-4804 + 14 x HCS-3804 | 7x2 + 7x2                          | 392            | 140              | 560   | 340/279 | 130                               | 133                      | 75                        |
| side opening   | 3 inch   | 14 x NIR-316U + 14 x NIR-312U | 7x2 + 7x2                          | 392            | 139              | 559   | 324/279 | 135                               | 135                      | 75                        |

## STAINLESS STEEL RACK SOLUTIONS

| Vertical rack type   | Box type | Rack/quantity<br>Stainless steel | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |        |        | Maximum box dimensions            |                          |                           |
|--|----------|----------------------------------|------------------------------------|----------------|------------------|--------|--------|-----------------------------------|--------------------------|---------------------------|
|  |          |                                  |                                    |                | Width            | Depth  | Height | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| <b>MDF-DU302VX-PE</b>  |          |                                  |                                    |                |                  |        |        |                                   |                          |                           |
| with trays   | 2 inch   | 12 x SDR-524-N                   | 3x4                                | 240            | 139.45           | 565.4  | 288.79 | 134                               | 138                      | 54                        |
| side opening   | 2 inch   | 12 x SUR-524-N                   | 3x4                                | 240            | 139.7            | 569.72 | 279.65 | 137                               | 138                      | 54                        |
| with trays   | 3 inch   | 12 x SDR-334-N                   | 3x4                                | 144            | 139.45           | 565.4  | 250.69 | 134                               | 138                      | 78                        |
| side opening   | 3 inch   | 12 x SUR-334-N                   | 3x4                                | 144            | 139.7            | 569.72 | 244.85 | 137                               | 138                      | 78                        |
| <b>MDF-DU502VX-PE / MDF-DU502VH-PE / MDF-DU502VH-PE / MDF-DU503VH-PE</b> |          |                                  |                                    |                |                  |        |        |                                   |                          |                           |
| with trays   | 2 inch   | 16 x SDR-624-N                   | 4x4                                | 384            | 139.45           | 565.4  | 325.12 | 134                               | 137                      | 52                        |
|  |          | 16 x SDR-624-P                   | 4x4                                | 384            | 139.45           | 565.4  | 339.59 | 134                               | 137                      | 54,5                      |
| side opening   | 2 inch   | 16 x SUR-624-N                   | 4x4                                | 384            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 52                        |
|  |          | 16 x SUR-624-P                   | 4x4                                | 384            | 139.7            | 569.72 | 332.74 | 136                               | 137                      | 54,5                      |
| with trays   | 3 inch   | 16 x SDR-434-N                   | 4x4                                | 256            | 139.45           | 565.4  | 320.54 | 134                               | 137                      | 78                        |
| side opening   | 3 inch   | 16 x SUR-434-N                   | 4x4                                | 256            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 78                        |
| <b>MDF-DU702VX-PE / MDF-DU702VH-PE / MDF-DU702VH-PE / MDF-DU703VH-PE</b> |          |                                  |                                    |                |                  |        |        |                                   |                          |                           |
| with trays   | 2 inch   | 24 x SDR-624-N                   | 6x4                                | 576            | 139.45           | 565.4  | 325.12 | 134                               | 137                      | 52                        |
|  |          | 24 x SDR-624-P                   | 6x4                                | 576            | 139.45           | 565.4  | 339.59 | 134                               | 137                      | 54,5                      |
| side opening   | 2 inch   | 24 x SUR-624-N                   | 6x4                                | 576            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 52                        |
|  |          | 24 x SUR-624-P                   | 6x4                                | 576            | 139.7            | 569.72 | 332.74 | 136                               | 137                      | 54,5                      |
| with trays   | 3 inch   | 24 x SDR-434-N                   | 6x4                                | 384            | 139.45           | 565.4  | 320.54 | 134                               | 137                      | 78                        |
| side opening   | 3 inch   | 24 x SUR-434-N                   | 6x4                                | 384            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 78                        |
| <b>MDF-DU901VHL-PE</b>   |          |                                  |                                    |                |                  |        |        |                                   |                          |                           |
| with trays   | 2 inch   | 28 x SDR-624-N                   | 7x4                                | 672            | 139.45           | 565.4  | 325.12 | 134                               | 137                      | 52                        |
| side opening   | 2 inch   | 28 x SUR-624-N                   | 7x4                                | 384            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 52                        |
| with trays   | 3 inch   | 14 x SDR-334-N                   | 7x2                                | 392            | 139.45           | 565.4  | 250.69 | 136                               | 137                      | 78                        |
|  |          | 14 x SDR-434-N                   | 7x2                                | 392            | 139.45           | 565.4  | 320.29 | 136                               | 137                      | 78                        |
| side opening   | 3 inch   | 14 x SUR-334-N                   | 7x2                                | 392            | 139.7            | 569.72 | 244.85 | 136                               | 137                      | 78                        |
|  |          | 14 x SUR-434-N                   | 7x2                                | 392            | 139.7            | 569.72 | 320.29 | 136                               | 137                      | 78                        |

\* Unit: mm

## TWINGUARD CHEST FREEZER INVENTORY RACKS

Models: MDF-DC500VX-PE | MDF-DC700VX-PE

### ALUMINIUM RACK SOLUTIONS

| Vertical rack type | Box type | Rack/quantity<br>Aluminium | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |        | Maximum box dimensions            |                          |                           |
|--------------------|----------|----------------------------|------------------------------------|----------------|------------------|-------|--------|-----------------------------------|--------------------------|---------------------------|
|                    |          |                            |                                    |                | Width            | Depth | Height | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| MDF-DC500VX-PE     |          |                            |                                    |                |                  |       |        |                                   |                          |                           |
| side opening       | 2 inch   | 32 x NIR-213C              | 8x4                                | 416            | 139              | 137   | 745    | 133                               | 133                      | 52                        |
| side opening       | 3 inch   | 32 x NIR-309C              | 8x4                                | 288            | 139              | 137   | 745    | 133                               | 133                      | 75                        |
| MDF-DC700VX-PE     |          |                            |                                    |                |                  |       |        |                                   |                          |                           |
| side opening       | 2 inch   | 40 x NIR-213C              | 10x4                               | 520            | 139              | 137   | 745    | 133                               | 133                      | 53                        |
| side opening       | 3 inch   | 40 x NIR-309C              | 10x4                               | 360            | 139              | 137   | 745    | 133                               | 133                      | 75                        |

### STAINLESS STEEL RACK SOLUTIONS

| Vertical rack type | Box type | Rack/quantity<br>Stainless steel | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |        | Maximum box dimensions            |                          |                           |
|--------------------|----------|----------------------------------|------------------------------------|----------------|------------------|-------|--------|-----------------------------------|--------------------------|---------------------------|
|                    |          |                                  |                                    |                | Width            | Depth | Height | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| MDF-DC500VX-PE     |          |                                  |                                    |                |                  |       |        |                                   |                          |                           |
| side opening       | 2 inch   | 32 x SCR-132-N                   | 8x4                                | 416            | 139.7            | 144   | 725.93 | 136                               | 142                      | 54                        |
| side opening       | 3 inch   | 32 x SCR-093-N                   | 8x4                                | 288            | 139.7            | 144   | 738.63 | 136                               | 142                      | 75                        |
| MDF-DC700VX-PE     |          |                                  |                                    |                |                  |       |        |                                   |                          |                           |
| side opening       | 2 inch   | 40 x SCR-132-N                   | 10x4                               | 520            | 139.7            | 144   | 725.93 | 136                               | 142                      | 54                        |
| side opening       | 3 inch   | 40 x SCR-093-N                   | 10x4                               | 360            | 139.7            | 144   | 738.63 | 136                               | 142                      | 75                        |

\* Unit: mm

## FROSTLESS FREEZERS INVENTORY RACKS

Models: MDF-DU500ZH-PE | MDF-DU700ZH-PE

### ALUMINIUM RACK SOLUTIONS

| Vertical rack type | Box type | Rack/quantity<br>Aluminium      | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |         | Maximum box dimensions            |                          |                           |
|--------------------|----------|---------------------------------|------------------------------------|----------------|------------------|-------|---------|-----------------------------------|--------------------------|---------------------------|
|                    |          |                                 |                                    |                | Width            | Depth | Height  | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| MDF-DU500ZH-PE     |          |                                 |                                    |                |                  |       |         |                                   |                          |                           |
| with trays         | 2 inch   | 8x HCS-5584SB + 8x HCS-6564SB   | 4x2 + 4x2                          | 352            | 139              | 560   | 290/339 | 133                               | 133                      | 53                        |
| side opening       | 2 inch   | 8x NIR-220U + 8x NIR-224U       | 4x2 + 4x2                          | 352            | 140              | 560   | 290/339 | 133                               | 133                      | 53                        |
| with trays         | 3 inch   | 8x HCS-3964SB + 8x HCS-4854SB   | 4x2 + 4x2                          | 224            | 139              | 560   | 290/339 | 133                               | 133                      | 90/80                     |
| side opening       | 3 inch   | 8x NIR-312U X 8x NIR-316U       | 4x2 + 4x2                          | 224            | 140              | 560   | 290/339 | 133                               | 133                      | 90/80                     |
| with trays         | 2 inch   | 4x HCS-519                      | 2x2                                | 352            | 280              | 560   | 640     | 133                               | 133                      | 53                        |
| MDF-DU700ZH-PE     |          |                                 |                                    |                |                  |       |         |                                   |                          |                           |
| with trays         | 2 inch   | 12x HCS-5584SB + 12x HCS-6564SB | 6x2 + 6x2                          | 528            | 139              | 560   | 290/339 | 133                               | 133                      | 53                        |
| side opening       | 2 inch   | 12x NIR-220U + 12x NIR-224U     | 6x2 + 6x2                          | 528            | 140              | 560   | 290/339 | 133                               | 133                      | 53                        |
| with trays         | 3 inch   | 12x HCS-3964SB + 12x HCS-4854SB | 6x2 + 6x2                          | 336            | 139              | 560   | 290/339 | 133                               | 133                      | 90/80                     |
| side opening       | 3 inch   | 12x NIR-312U X 12x NIR-316U     | 6x2 + 6x2                          | 336            | 140              | 560   | 290/339 | 133                               | 133                      | 90/80                     |
| with trays         | 2 inch   | 6x HCS-519                      | 3x2                                | 528            | 280              | 560   | 640     | 133                               | 133                      | 53                        |

### STAINLESS STEEL RACK SOLUTIONS

| Vertical rack type | Box type | Rack/quantity<br>Aluminium    | Freezer layout<br>(Columns x Rows) | Total<br>boxes | Rack dimensions* |       |         | Maximum box dimensions            |                          |                           |
|--------------------|----------|-------------------------------|------------------------------------|----------------|------------------|-------|---------|-----------------------------------|--------------------------|---------------------------|
|                    |          |                               |                                    |                | Width            | Depth | Height  | Footprint<br>of box<br>max width* | Lid of box<br>max width* | Box height<br>max height* |
| MDF-DU500ZH-PE     |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays         | 2 inch   | 8x SDR-624-P + 8x SDR-524-P   | 4x2 + 4x2                          | 352            | 139              | 565   | 325/289 | 134                               | 137                      | 52/54                     |
| side opening       | 2 inch   | 8x SUR-624-P + 8x SUR-524-P   | 4x2 + 4x2                          | 352            | 140              | 570   | 320/280 | 136                               | 137                      | 52/54                     |
| with trays         | 3 inch   | 8x SDR-434-N + 8x SDR-334-N   | 4x2 + 4x2                          | 224            | 139              | 565   | 321/251 | 134                               | 137                      | 78                        |
| side opening       | 3 inch   | 8x SUR-434-N + 8x SUR-334-N   | 4x2 + 4x2                          | 224            | 140              | 570   | 320/245 | 136                               | 137                      | 78                        |
| MDF-DU700ZH-PE     |          |                               |                                    |                |                  |       |         |                                   |                          |                           |
| with trays         | 2 inch   | 12x SDR-624-P + 12x SDR-524-P | 6x2 + 6x2                          | 528            | 139              | 565   | 325/289 | 134                               | 137                      | 52/54                     |
| side opening       | 2 inch   | 12x SUR-624-P + 12x SUR-524-P | 6x2 + 6x2                          | 528            | 140              | 570   | 320/280 | 136                               | 137                      | 52/54                     |
| with trays         | 3 inch   | 12x SDR-434-N + 12x SDR-334-N | 6x2 + 6x2                          | 336            | 139              | 565   | 321/251 | 134                               | 137                      | 78                        |
| side opening       | 3 inch   | 12x SUR-434-N + 12x SUR-334-N | 6x2 + 6x2                          | 336            | 140              | 570   | 320/245 | 136                               | 137                      | 78                        |

\* Unit: mm

# PRO ECO / VIP ECO INVENTORY RACKS

Models: MDF-DU300H-PE | MDF-DC102VH-PE | MDF-DC202VH-PE

## ALUMINIUM RACK SOLUTIONS

| Vertical rack type         | Box type | Rack/quantity Aluminium                      | Freezer layout (Columns x Rows) | Total boxes | Rack dimensions* |         |         | Maximum box dimensions      |                       |                        |
|----------------------------|----------|--|---------------------------------|-------------|------------------|---------|---------|-----------------------------|-----------------------|------------------------|
|                            |          |  |                                 |             | Width            | Depth   | Height  | Footprint of box max width* | Lid of box max width* | Box height max height* |
| <b>MDF-DU300H-PE</b>       |          |  |                                 |             |                  |         |         |                             |                       |                        |
| with trays                 | 2 inch   | 6 x HCS-32-4584/143<br>+ 6 x HCS-32-5584/143 | 3x2 + 3x2                       | 216         | 143              | 560     | 232/280 | 133                         | 136                   | 52                     |
| side opening               | 2 inch   | 6 x NIR-216U<br>+ 6 x NIR-220U               | 3x2 + 3x2                       | 216         | 139              | 559     | 232/279 | 135                         | 135                   | 52                     |
| with trays                 | 3 inch   | 12 x HCS-32-3804/143                         | 3x4                             | 144         | 143              | 560     | 232     | 133                         | 136                   | 71                     |
| side opening               | 3 inch   | 12 x NIR-312U                                | 3x4                             |             | 144              | 139     | 559     | 279                         | 135                   | 135                    |
| <b>MDF-DC102VH-PE</b>      |          |  |                                 |             |                  |         |         |                             |                       |                        |
| side opening               | 2 inch   | 6 x NIR-207C                                 | 3x2                             | 42          | 142              | 141     | 403     | 133                         | 133                   | 53                     |
| ladder type                | 2 inch   | 21 x VC-660/3ST                              | 6x3 + 3                         | 63          | 58               | 133     | 424     | 133                         | 133                   | 53                     |
| side opening               | 3 inch   | 6 x NIR-305C                                 | 3x2                             | 30          | 142              | 141     | 403     | 133                         | 133                   | 75                     |
| side opening + ladder type | 2 inch   | 6 x NIR-207C + 5x VC-660/3ST                 | 3x2 + 5                         | 57          | 142/58           | 141/133 | 403/424 | 133                         | 133                   | 53                     |
| <b>MDF-DC202VH-PE</b>      |          |  |                                 |             |                  |         |         |                             |                       |                        |
| side opening               | 2 inch   | 18 x NIR-207C                                | 6x3                             | 126         | 142              | 141     | 403     | 133                         | 133                   | 53                     |
| ladder type                | 2 inch   | 47 x VC-660/3ST                              | 13x3 + 8                        | 141         | 58               | 133     | 424     | 133                         | 133                   | 53                     |
| side opening               | 3 inch   | 18 x NIR-305C                                | 6x3                             | 90          | 142              | 141     | 403     | 133                         | 133                   | 75                     |

## STAINLESS STEEL RACK SOLUTIONS

| Vertical rack type    | Box type | Rack/quantity Stainless steel    | Freezer layout (Columns x Rows) | Total boxes | Rack dimensions* |        |                 | Maximum box dimensions      |                       |                        |
|-----------------------|----------|----------------------------------|---------------------------------|-------------|------------------|--------|-----------------|-----------------------------|-----------------------|------------------------|
|                       |          |                                  |                                 |             | Width            | Depth  | Height          | Footprint of box max width* | Lid of box max width* | Box height max height* |
| <b>MDF-DU300H-PE</b>  |          |                                  |                                 |             |                  |        |                 |                             |                       |                        |
| with trays            | 2 inch   | 6 x SDR-424-N<br>+ 6 x SDR-524-N | 3x2 + 3x2                       | 216         | 139.45           | 565.4  | 288.79 / 325.12 | 134                         | 137                   | 54                     |
| side opening          | 2 inch   | 6 x SUR-424-N<br>+ 6 x SUR-524-N | 3x2 + 3x2                       | 216         | 139.7            | 569.72 | 233.68 / 279.65 | 134                         | 137                   | 54                     |
| with trays            | 3 inch   | 12 x SDR-334-N                   | 3x4                             | 144         | 139.45           | 565.4  | 250.69          | 134                         | 137                   | 78                     |
| side opening          | 3 inch   | 12 x SUR-334-N                   | 3x4                             | 144         | 139.7            | 569.72 | 244.85          | 136                         | 137                   | 78                     |
| <b>MDF-DC102VH-PE</b> |          |                                  |                                 |             |                  |        |                 |                             |                       |                        |
| side opening          | 2 inch   | 6 x SCR-072-N                    | 3x2                             | 42          | 139.7            | 144    | 397             | 136                         | 142                   | 54                     |
| side opening          | 3 inch   | 6 x SCR-053-N                    | 3x2                             | 30          | 139.70           | 144    | 414.78          | 136                         | 142                   | 75                     |
| <b>MDF-DC202VH-PE</b> |          |                                  |                                 |             |                  |        |                 |                             |                       |                        |
| side opening          | 2 inch   | 18 x SCR-072-N                   | 6x3                             | 126         | 139.7            | 144    | 397             | 136                         | 142                   | 54                     |
| side opening          | 3 inch   | 18 x SCR-053-N                   | 6x3                             | 90          | 139.70           | 144    | 414.78          | 136                         | 142                   | 75                     |

\* Unit: mm







# VALIDATION & QUALIFICATION SOLUTIONS

PHC Europe BV is a vertical component manufacturer that can provide turn-key solutions for validation and qualification in accordance with all current GMPs, GLPs, GCPs, 21 CFR Part 11, PAT, ISO and specific customer requirements and applications. Because many of our key component parts are designed and built by PHC Europe BV, we offer the most precise and in-depth validation resources specific to PHCbi laboratory products. Whatever your validation needs are, PHCbi provides comprehensive expertise in laboratory equipment to meet your exact compliance needs. PHCbi validation systems employ advanced technology coupled with the latest trends to insure compliance with accurate and time efficient completion.

## Validation & Qualification Solutions for laboratory equipment

Turn key solutions available for:

- Ultra-Low Freezers
- Cryogenic Freezers
- Biomedical Freezers
- Bloodbank refrigerators
- Pharmaceutical refrigerators
- Incubators
- Ovens
- Autoclaves
- Environmental test chambers



# Installation and Operational Qualification

## Qualification IOQ

PHC Europe BV offer onsite validation of PHCBI supplied equipment via Installation and Operational Qualification Protocol IOQ.

## Installation Qualification (IQ)

Verifies and documents the equipment installation is compliant with the manufacturer's requirements and specifications.

## Operational Qualification (OQ)

Verifies and documents the full functional operation of the installed equipment (as specified by PHCBI or other OEM supplied equipment). Temperature performance will be mapped over a continuous 24hr period; also, a short open-door test included toward the end of this period. Data produced will be compared with manufacturers published equipment specification. Product specific parameters such as CO<sub>2</sub>/O<sub>2</sub>, %RH etc. are included within the relevant equipment IOQ protocol.

## Additional options:


Additional temperature sensor positions, Extended logging period 48/72hrs, Simulated customer loaded mapping, Powerfail/ Recovery performance tests.

**Process Qualification (PQ)** is usually conducted and performed by customer/end-user as the equipment will be in an environment where specific user conditions and actual product is loaded, stored and accessed e.g. the customers actual production or product processing area. The PQ will probably refer to customer specific Standard Operating Procedures documents (SOPs).

| PHCbi validation services by model                               |       |                 |                |     |                    |          |
|--|-------|-----------------|----------------|-----|--------------------|----------|
|  | Temp. | CO <sub>2</sub> | O <sub>2</sub> | %RH | Lighting (Lux/Par) | Pressure |
| Validation by Model  | ✓     |                 |                |     |                    |          |
| MDF-150°C Freezers   | ✓     |                 |                |     |                    |          |
| MDF-86°C Freezers  | ✓     |                 |                |     |                    |          |
| MDF-30°C Freezers  | ✓     |                 |                |     |                    |          |
| MBR Blood Bank Refrigerators                                     | ✓     |                 |                |     |                    |          |
| MPR Pharmaceutical Refrigerators                                 | ✓     |                 |                |     |                    |          |
| MIR Incubator Series   | ✓     |                 |                |     |                    |          |
| MCO CO <sub>2</sub> & O <sub>2</sub> /CO <sub>2</sub> Incubators | ✓     | ✓               | ✓              | ✓   |                    |          |
| MLS Top Loading Autoclaves                                       | ✓     |                 |                |     |                    | ✓        |
| MLR Environmental Test Chamber                                   | ✓     |                 |                | ✓   | ✓                  |          |
| CBS Standard LN <sub>2</sub> Freezers                            | ✓     |                 |                |     |                    |          |

Example: product identification and specific storage requirements; loading patterns etc. therefore making the PQ a unique and customer specific document. PHC Europe however will provide assistance to customers where required in either the preparation or assisted execution of the Process Qualification.

PHC Europe BV are also able to offer a "Temperature Mapping Service" for customers wishing to verify actual equipment performance as installed, this service is also available for all NON PHCBI equipment.





# PHC Europe

A Member of PHC Group

Eikdonk 1 | 4825 AZ Breda | Netherlands  
T: +31 (0) 76 543 3833

[www.phcd.com/eu/biomedical](http://www.phcd.com/eu/biomedical)

**PHCbi**